# Respirator PFE decision rule guide

This document is a guide on PFE testing criteria for analysts and validators. It is imperative that the correct PFE class is used, matching the highest claimed PFE claimed by the manufacture when using this decision-making guide. This document should be used alongside the SOP "Particulate Penetration Efficiency Assessment – Rapid Screening" found <u>D20-3958656</u>.

### Initial 10 specimens (specimens 1-10)

If all specimens tested pass – cease testing, sample passes.

If only 1 specimen is within the lower limit of acceptance, another 10 specimens of the same sample must be tested to determine the overall result – continue testing.

If 1 or more specimens less than the lower limit of acceptance – cease testing, sample fails.

If 2 or more specimens within the lower limit of acceptance – cease testing, sample fails.

## Follow-up 10 specimens (specimens 11-20)

Test another 10 specimens from the same sample, only if 1 specimen was within the lower limit of acceptance.

If all specimens tested pass – cease testing, sample passes.

If only 1 specimen within the lower limit of acceptance (2 in total) – cease testing, sample passes.

If 1 or more specimens less than the lower limit of acceptance – cease testing, sample fails.

If 2 or more specimens within the lower limit of acceptance (≥3 in total) – cease testing, sample fails.

#### Initial 10 specimens (specimens 1-10)

PFE Class	Pass	Indeterminate	Minor Fail	Major Fail
	first 10 specimens	1 specimen within	≥2 specimen	≥ 1 specimen
	above limit	lower limit of	within lower	<b>below</b> lower
		acceptance	limit of	limit of
			acceptance	acceptance
P1	≥ 80 %	79-80 %	79-80 %	< 79 %
KN90	≥ 90 %	89-90 %	89-90 %	< 89 %
P2/FFP2/KF94	≥ 94 %	93-94 %	93-94 %	< 93 %
KN95/N95	≥ 95 %	94-95 %	94-95 %	< 94 %
FFP3/N99	≥ 99 %	98-99 %	98-99 %	< 98 %
P3	≥ 99.95%	98.95-99.95%	98.95-99.95%	< 98.95 %

#### Follow up 10 specimens (specimens 11-20)

PFE Class	Pass next 10 specimens above limit	Pass 1 follow-up specimen, 2 total, within lower limit of acceptance	Minor Fail ≥3 total speicmens within lower limit of acceptance	Major Fail ≥ 1 specimen below lower limit of acceptance
P1	≥ 80 %	79-80 %	79-80 %	< 79 %
KN90	≥ 90 %	89-90 %	89-90 %	< 89 %
P2/FFP2/KF94	≥ 94 %	93-94 %	93-94 %	< 93 %
KN95/N95	≥ 95 %	94-95 %	94-95 %	< 94 %
FFP3/N99	≥ 99 %	98-99 %	98-99 %	< 98 %
P3	≥ 99.95%	98.95-99.95%	98.95-99.95%	< 98.95 %

## **Analyst decision table**

	Failures or lower limit of acceptance	Action	Failures or lower limit of acceptance	Result
Test Initial 10	0		<b>—</b>	Sample Passes
	≥ 1 Major Failure ≥ 2 Minor Failure			Sample Fails
	1 within lower limit of acceptance	Test another 10	1 lower limit of acceptance	Sample Passes (Total of 2 within lower limit of acceptance)
			≥ 1 Major Failure ≥ 2 Minor Failures (in this group of 10)	Sample Fails

## Analyst decision tree

