



## Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 31 - ID : 420887

Released by **S22** on 21/11/2018 10:36:59

Report #: <input type="text" value="56421"/>	Records Management #: <input type="text"/>	Reporter's Reference #: <input type="text" value="0703058556"/>	Report Type: <input type="text" value="Final"/>
ARTG: 293808	<a href="#">Document Container URL</a>		

Report Information Section

Report Status: <input type="text" value="Closed"/>	Sponsor's Reported Category: <input type="text" value="Death / Serious Injury"/>	Date of Adverse Event: <input type="text" value="S22"/>	Date of Initial Report: <input type="text" value="19/03/2019"/>
Date of Final Report: <input type="text" value="05/04/2019"/>	Date of Initial TGA Act on: <input type="text" value="19/03/2019"/>	Reviewed by Team: <input type="text"/>	Date Response Received: <input type="text"/>
Date Completed: <input type="text" value="08/04/2019"/>	Operator at Time of Event: <input type="text" value="Healthcare Professional"/>	If 'Other' Operator Selected: <input type="text"/>	Reporter Confidentiality: <input type="text" value="No"/>
Source of Report: <input type="text" value="Sponsor"/>	If 'Other' Source Selected: <input type="text"/>	Type of Initial Action: <input type="text" value="Trend data only"/>	

### Event Description for Website Publication:

**An infection from anastomotic leak was noted.**

### Clinical Event Information:

According to the reporter, four days postoperative from a laparoscopic high anterior resection, an infection from anastomotic leak was noted. It is unknown if the leak is related to the stapler. The clinical anastomotic leak was determined by CT scan and regal contrast leaked out of anastomosis. There were no issues during the operation. The stapler appeared to fire correctly, it passed an air leak test after the stapler had been used, and normal full tissue donuts were seen. The patient had initially undergone Hartmann's procedure and hospitalisation was extended for one week due to the leak.

Number of Incidents in Report: <input type="text" value="1"/>	Contact: <input type="text"/>	Alternative Person Title: <input type="text"/>	Alternative Person First Name: <input type="text"/>
Alternative Person Surname: <input type="text"/>	Alternative Person Phone: <input type="text"/>	Alternative Person Fax: <input type="text"/>	Alternative Person Email: <input type="text"/>

### Patient Information

Sex: <input type="text" value="S22"/>	Weight: <input type="text" value="Unknown"/>	Age: <input type="text" value="S22"/>
Patient Focused Corrective Action Taken: <input type="text" value="Hospitalisation was extended by one week due to the leak."/>	Patient History: <input type="text" value="Crohn's disease and fistula between colon and bladder."/>	
Patient Outcome/Consequences: <input type="text" value="Alive - with injury, requiring extended hospitalisation."/>	Additional Event Description: <input type="text"/>	
Describe any test (Lab, xray, etc.): <input type="text"/>	Injured - Extent of Injury: <input type="text" value="Temporary Injury"/>	Consequence: <input type="text"/>
Medical Problem Device Used For: <input type="text"/>	Additional Patients Added: <input type="text" value="0"/>	Other medical devices currently using/implanted: <input type="text"/>

### Submitting Reporter Section

Search Reporter By Surname: <input type="text" value="S22"/>	Reporter #: <input type="text"/>	Preferred Contact Method: <input type="text"/>
Reporter Title: <input type="text" value="S22"/>	First Name: <input type="text" value="S22"/>	Surname: <input type="text" value="S22"/>
Position: <input type="text" value="Post Market Vigilance Specialist"/>	Company/Institution: <input type="text" value="Medtronic Australasia Pty Ltd."/>	
Address 1: <input type="text"/>	Address 2: <input type="text"/>	Town/Suburb: <input type="text"/>
		State: <input type="text"/>

2 Alma Road  
 Country:  
 Australia  
 Mobile:

2 Alma Road  
 Postcode:  
 2113  
 Email:  
 s22@medtron c.com

Macquarie Park  
 Phone:  
 s22  
 Are you happy for the device company to contact you about the incident?:

NSW **Document 1**  
 Fax:  
 Last External Submission By:  
 106402\_837 - 05/04/2019 12:03

## Initial Reporter Section

As Above?:  
 No  
 Search Reporter By Surname:  
 s22  
 Title:  
 Position:  
 Address 1:  
 Ipswich Road  
 Postcode:  
 4102  
 Mobile:

If No, fill out the following:

Initial Reporter #:  
 First Name:  
 s22  
 Surname:  
 s22  
 Company/Inst tution:  
 PRINCESS ALEXANDRA HOSPITAL  
 Town/Suburb:  
 Brisbane  
 Phone:  
 Allow the dev ce company to contact you about the incident:

Initial Reporter Conf dential:  
 Yes  
 Preferred Contact Method:

## Dev ce Information Sect on

Product Exempt (Note: If not exempt, enter ARTG No):  
 No  
 Product Licence Category:  
 Included  
 Brand Name:  
 Covidien EEA Circular Stapler - Intraluminal circular stapler, single-use  
 Model #:  
 TRIEEA28MT  
 Purchase Date:  
 Place of Implantation:  
 Access Contact Surname:  
 Add tional Dev ces Added:  
 0

Search Device ARTG:  
 293808  
 Dev ce Class:  
 Class IIb  
 In tial Dev ce Descript on:  
 Covidien EEA Circular Stapler - Intraluminal circular stapler, single-use  
 Serial #:  
 Expiry Date:  
 Reported Device Location:  
 Discarded  
 Access Contact Phone:

Dev ce ARTG #:  
 293808  
 GMDN / UMDN Code:  
 59875  
 Usage of Device:  
 Single Use  
 Batch #:  
 Date of Implant:  
 Access Contact Title:  
 Access Contact Fax:

Therapeutic L cence Type:  
 Medical Dev ce  
 GMDN / UMDN Text:  
 Intraluminal circular stapler, single-use  
 Software Vers on:  
 Lot #:  
 Unknown  
 Date of Explant:  
 Access Contact First Name:  
 Access Contact Email:

## Manufacturer Informat on Section

Manufacturer Name:  
 Covidien llc  
 Address 2:  
 Postcode:  
 Town/Suburb:  
 Mansfield  
 Phone:

Manufacturer Client Id:  
 54968  
 State/Province:  
 MA  
 Fax:

Address 1:  
 15 Hampshire Street  
 Country:  
 United States  
 Email:

02048			
Manufacturer Informed:	Date Aware of Adverse Event:	Contact Title:	Contact First Name:
Yes	15/03/2019		
Contact Surname:			

Supplier Information Section

Supplier Name:	Address 1:		Address 2:
Town/Suburb:	State:	Country:	Postcode:
Phone:	Fax:	Email:	Website:
Supplier Informed:	Date of Supplier Contact:	Contact Title:	Contact First Name:
Contact Surname:	Contact Phone:	Contact Fax:	Contact Email:

Report Information - duplicated information from other parts of the report, for use in risk assessments.

Licence Start Date:	Date of Initial TGA Action:	Report Status:	
13/09/2017	19/03/2019	Closed	
Problems Observed:			
; ; Mechanical Problem; Leak / Splash;			

Report Status				
For website publication:	Ready for Publication:	Investigated:	Investigation Reason:	Team Assignment:
Yes	Yes	No	Device not returned	Team B (Iib, Is & Im)
Report Priority:				
Not Investigated				

Team Review				
Reviewed by Team:	Reason Sent To Meeting:	Outcome from team meeting:		
Team Meeting Notes:				

DPRC Review				
Reviewed by DPRC:	DPRC Reason Sent To Meeting:	Outcome from DPRC Meeting:		
Meeting Notes:				

Initial Risk Analysis				
Date:	Assessor:	Licence Status:	Status Reason:	Status Effective Date:
20/03/2019	S22	Active		13/09/2017
Injured Party:	Potential Effect:	Actual Effect:	Found Prior To Use:	Sample Received:
Patient	Serious Injury	Temporary Injury	No	No
Sterile:	Invasive Device:	Single Use:	Human Origin:	Genetically Modified:
Yes	Yes	Yes	No	No
Reusable:	Risk Frequency:	Risk Severity:	Risk Rating:	Further Review Needed:

No  Serious

Risk Assessment Notes:

Final Risk Assessment:

No

RISK RATING	Severity				
	Life-threatening	Serious	Minor	Nil	Unknown
Frequently	Critical Risk	Critical Risk	Major Risk	Minor Risk	Major Risk
Sometimes	Critical Risk	Major Risk	Minor Risk	Minor Risk	Minor Risk
Rarely	Major Risk	Minor Risk	Minor Risk	Non-significant Risk	Minor Risk
Unlikely	Minor Risk	Minor Risk	Non-significant Risk	Non-significant Risk	Non-significant Risk
Unknown	Major Risk	Minor Risk	Minor Risk	Non-significant Risk	No risk assessment

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section				
Update Device Details?:	Date:	Assessor:	License Status:	Status Reason:	Status Effective Date:
Yes	08/04/2019	§22	Active	Active	13/09/2017
Copy Data From:	Injured Party:	Potential Effect:	Actual Effect:	Found Prior To Use:	Sample Received:
In trial	Patient	Serious Injury	Temporary Injury	No	No
	Sterile:	Invasive Device:	Single Use:	Human Origin:	Genetically Modified:
	Yes	Yes	Yes	No	No
	Reusable:	Risk Frequency:	Risk Severity:	Risk Rating:	Final Risk Assessment:
	No	Rarely	Serious	Minor Risk	Yes
	Risk Assessment Notes:				
	<input type="text"/>				

Sponsor/Manufacturer Information Section

Search Sponsors:	Name:	Client #:
837	Medtronic Australasia Pty Ltd	837
Attention To:	Address 1:	Town/Suburb:
§22	PO Box 945	NORTH RYDE BC
State:	Postcode:	Phone:
NSW	1670	§22
Email:		Fax:
§22@medtronic.com		

Investigation Information Section - Submitted by Sponsor/Manufacturer

Device Analysis Results:	Corrective/Preventative Actions:
The reported EEA* Circular Stapler with Tri-Staple* Technology 28mm Medium/Thick sample or additional supporting materials from the account was not available for this incident. Therefore, a definitive root cause could not be determined with regard to the reported condition. Records from each manufacturing lot are thoroughly reviewed to ensure that products are released meeting all Medtronic quality release specifications at the time of manufacture. Should new information become available, the file will be re-opened and the investigation will be amended as appropriate.	There are no actions deemed necessary at this time.
Details of Similar Events:	Additional Details (use for tables):
<input type="text"/>	<input type="text"/>

Nil other similar incidents in Australia.  
 Similar incidents worldwide: 2/1,325. Austria 1, Germany 1

CAPA# Reference:



Risk Assessment

Frequency:  Severity:

Rating:

Type Cause and Outcome:  Number of Similar Events:

Expected Rate:  Actual Rate:

Countries Similar Events Also Occurred:  
 Austria 1, Germany 1

Completed Actions:  Planned Actions and Proposed Timelines:

Additional Comments:

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes
<input type="checkbox"/>	Final Report Due	TGA Regulatory Action	Regulatory Information Request (s41)		19/03/2019	27/06/2019			

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Mechanical Problem	Leak / Splash		Infection

Investigation Findings

Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
No Findings Available			

Investigation Conclusion

Investigation Conclusion (L1)	Investigation Conclusion (L2)	If Additional Conclusion Detail Requested
Cause Not Established		

Investigation Outcomes

Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested
Reviewed, for Trending Purposes Only		

Investigation Summary

Investigation Type: <input type="text"/>	Latest Investigat on (DII) where this DIR is the Primary DIR: <input type="text"/>	Latest Investigat on (DII) where this DIR is a Related DIR: <input type="text"/>	Investigator: <input type="text"/>	Extens on Number: <input type="text"/>
Investigator's Notes: <input type="text"/>		Summary Findings: No further investigation will occur at this time, however the TGA will continue to mon tor the rate and pattern of occurrence and may re-open the file as appropriate.		Recall Number: <input type="text"/>

**Note:** Letter generat on buttons disabled if report not ready for webs te publ cation or risk analysis not completed.

Device Lookup

This section is used to match informat on prov ded via UDIR forms to ARTG information. You can select a Brand/Name from informat on prov ded in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered): <input type="text"/>	Brand Name: <input type="text"/>	Manufacturer Name: <input type="text"/>	Device ARTG #: <input type="text"/>
---	-------------------------------------	--	--

Other Devices

Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Cl ck **New** to begin entering informat on.

Rec No	
1	

Samples Record - Click **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details			Add t onal Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Dev ce Descript on:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Vers on Number:	
						Who sent the dev ce to the TGA?:		Why does the TGA have the sample?:	

Additional Patients

Cl ck **[N]** to begin entering information.

Patient Details			
Sex:	Weight:	Age:	
Patient Focused Corrective Act on Taken:		Patient History:	

Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence: <b>Document 1</b>
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

Additional Device Information

Where did you get this device from?:  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

Similar Events

Similar events - how many times?:  Date of Recent Report:  Event Reported To:  Reporter Reference Number:

Device Access - Alternate Device Contact Information Provided

Title:  First Name:  Last Name:  Phone:

Fax:  Email:

Incident Location Details

Occurred in Australia:  Organisation:  Address Line 1:  Address Line 2:

Town/Suburb:  State:  Postcode:

Flow Details DIR-REQ - Device Incident Request 165127

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
165127	DIR-REQ		Closed	<b>S22</b>	OPR Administration User	08/04/2019	Normal	0

Signature Details

Role	IRIS Investigator
User	<b>S22</b>
Signed At	08/04/2019 08:44:22
Comment	



## Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 36 - ID : 459466

26/11/2019

SIGNED

Released by Theta Technologies on 13/02/2020 10:09:06

Report #:	Records Management #:	Reporter's Reference #:	Report Type:
<input type="text" value="60703"/>	<input type="text"/>	<input type="text" value="0703474545"/>	<input type="text" value="Final"/>
ARTG: 293808	<a href="#">Document Container URL</a>		

### Report Information Section

Report Status:	Sponsor's Reported Category:	Date of Adverse Event:	Date of Initial Report:
<input type="text" value="Closed"/>	<input type="text" value="Death / Serious Injury"/>	<input type="text" value="§22"/>	<input type="text" value="26/11/2019"/>
Date of Final Report:	Date of Initial TGA Action:	Reviewed by Team:	Date Response Received:
<input type="text" value="24/02/2020"/>	<input type="text" value="26/11/2019"/>	<input type="text" value="03/03/2020"/>	<input type="text" value="27/11/2019"/>
Date Completed:	Operator at Time of Event:	If 'Other' Operator Selected:	Reporter consents to contact by sponsor:
<input type="text" value="04/03/2020"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="N/A"/>
Source of Report:	If 'Other' Source Selected:	Type of Initial Action:	
<input type="text" value="Sponsor"/>	<input type="text"/>	<input type="text" value="For IRIS Meeting"/>	

### Event Description for Website Publication:

### Clinical Event Information:

According to the reporter, post-operatively on a robotic laparoscopic anterior resection for endometriosis procedure, after the stapler was fired and air leak test was performed and the donuts were inspected. it was also reported that the surgical procedure was extended for more than 30 minutes. It was reported that there were no air leaked and the donuts looked good. 4 days after the surgery the anastomosis dehiscence. There was a tissue loss and tissue damaged due to the reported event. Patient has undergone a procedure to create diverting ileostomy and another to repair a pelvic abscess that eroded into the bladder and vagina.

Number of Incidents in Report:	Contact:	Alternative Person Title:	Alternative Person First Name:
<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Alternative Person Surname:	Alternative Person Phone:	Alternative Person Fax:	Alternative Person Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

### Recorded Problems Observed

#### Recorded Problems Observed:

### Clinical Signs, Symptoms and Conditions

#### Recorded Clinical Signs, Symptoms and Conditions:

### Health Impact

#### Recorded Health Impacts:

### Patient Information

Sex:	Weight:	Age:
<input type="text" value="§22"/>	<input type="text"/>	<input type="text"/>
Patient Focused Corrective Action Taken:	Patient History:	
<input type="text" value="The patient spent 25 days in the hospital before being released. The patient underwent a Hartmann's reversal procedure."/>	<input type="text" value="Non-Smoker."/>	
Patient Outcome/Consequences:	Additional Event Description:	
<input type="text" value="Alive."/>	<input type="text"/>	
Describe any test (Lab, xray, etc.):	Injured - Extent of Injury:	Consequence:
<input type="text"/>	<input type="text"/>	<input type="text"/>
		Other medical devices currently using/implanted:
		<input type="text"/>



<input type="text"/>	<input type="text" value="Serious Injury"/>	<input type="text" value="Required surgical intervention"/>	<input type="text"/>
Medical Problem Device Used For:	Additional Patients Added:		
<input type="text"/>	0		

## Submitting Reporter Section

Search Reporter By Surname:	Reporter #:	Preferred Contact Method:	
<input type="text" value="§22"/>	<input type="text"/>	<input type="text"/>	
Reporter Title:	First Name:	Surname:	
<input type="text" value="§22"/>	<input type="text" value="§22"/>	<input type="text" value="§22"/>	
Position:	Company/Institution:		
<input type="text" value="§22"/>	<input type="text" value="Medtronic Australasia"/>		
Address 1:	Address 2:	Town/Suburb:	State:
<input type="text" value="2 Alma Road"/>	<input type="text"/>	<input type="text" value="Macquarie Park"/>	<input type="text" value="NSW"/>
Country:	Postcode:	Phone:	Fax:
<input type="text" value="Australia"/>	<input type="text" value="2113"/>	<input type="text" value="§22"/>	<input type="text"/>
Mobile:	Email:	Last External Submission By:	
<input type="text"/>	<input type="text" value="§22 @medtronic.com"/>	<input type="text" value="§22 - 24/02/2020 12:54"/>	

## Initial Reporter Section

As Above?:	<i>If No, fill out the following:</i>	Initial Reporter Confidential:	
<input type="text" value="No"/>		<input type="text" value="Yes"/>	
Search Reporter By Surname:	Initial Reporter #:	Preferred Contact Method:	
<input type="text" value="§22"/>	<input type="text"/>	<input type="text"/>	
Title:	First Name:	Surname:	
<input type="text"/>	<input type="text" value="§22"/>	<input type="text" value="§22"/>	
Position:	Company/Institution:		
<input type="text"/>	<input type="text" value="SYDNEY ADVENTIST HOSPITAL"/>		
Address 1:	Address 2:	Town/Suburb:	State:
<input type="text" value="Theatres, Lvl 5"/>	<input type="text"/>	<input type="text" value="WAHROONGA"/>	<input type="text" value="NSW"/>
Postcode:	Country:	Phone:	Fax:
<input type="text" value="2076"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Mobile:	Email:	Allow the device company to contact you about the incident:	
<input type="text"/>	<input type="text"/>	<input type="checkbox"/>	

## Device Information Section

Product Exempt ( <i>Note: If not exempt, enter ARTG No</i> ):	Search Device ARTG:	Device ARTG #:	Therapeutic Licence Type:
<input type="text" value="No"/>	<input type="text" value="293808"/>	<input type="text" value="293808"/>	<input type="text" value="Medical Device"/>
Product Licence Category:	Device Class:	GMDN / UMDN Code:	GMDN / UMDN Text:
<input type="text" value="Included"/>	<input type="text" value="Class IIb"/>	<input type="text" value="59875"/>	<input type="text" value="Intraluminal circular stapler, single-use"/>
Brand Name:	Initial Device Description:	Usage of Device:	Software Version:
<input type="text" value="EEA - Intraluminal circular stapler, single-use"/>	<input type="text" value="EEA - Intraluminal circular stapler, single-use"/>	<input type="text"/>	<input type="text"/>
Model #:	Serial #:	Batch #:	Lot #:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Purchase Date:	Expiry Date:	Date of Implant:	Date of Explant:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Date of Inital Procedure:	Place of Implantation:	Reported Device Location:	Access Contact Title:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Access Contact First Name:	Access Contact Surname:	Access Contact Phone:	Access Contact Fax:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Access Contact Email:	Licence Status:	Status Effective Date:	Additional Devices Added:
<input type="text"/>	A	13/09/2017	0

Manufacturer Information Section

Manufacturer Name:	Manufacturer Client Id:	Address 1:
Covidien llc	54968	<input type="text"/>
Address 2:	State/Province:	Country:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Postcode:	Fax:	Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Manufacturer Informed:	Contact Title:	Contact First Name:
Yes	<input type="text"/>	<input type="text"/>
Contact Surname:		
<input type="text"/>		

Supplier Information Section

Supplier Name:	Address 1:	Address 2:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Town/Suburb:	Country:	Postcode:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Phone:	Email:	Website:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Supplier Informed:	Contact Title:	Contact First Name:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Contact Surname:	Contact Fax:	Contact Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>

Report Status

For website publication:	Ready for Publication:	Investigated:	Investigation Reason:	Team Assignment:	Team Priority:
Yes	Yes	No	Device not returned	Team B (Iib, Is & Im)	Not Investigated

Team Review

Reviewed by Team:	Reason Sent To Meeting:	Outcome from team meeting:
03/03/2020	Not enough information has been provided	No further action
Notes for Team meeting:		
send to meeting due to level 1 investigation - device not returned. no similar DIR in 6 month. no recall/pmr		
Outcomes from Team Meeting:		
NFA - Device not returned. No similar DIRs.		

DPRC Review

Reviewed by DPRC:	DPRC Reason Sent To Meeting:	Outcome from DPRC Meeting:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Meeting Notes:		
<input type="text"/>		

Initial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Date: 27/11/2019	Severity: 5 - An illness/injury was resolved or prevented with treatment by a health professional	Incidents in the last 12 months:	Manufacturer analysis:	Assessor: §22	Manufacturer documentation: Unknown - updated information from the manufacturer is required
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION: Level 1 Investigation (to complete screening)	FINAL LEVEL OF INVESTIGATION: Screening only	Injured Party: Patient	Device Recalls: 0. No recalls for similar incidents in Australia
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL: Awaiting final report		Found Prior To Use: No	Is AE covered by current recall:
Incidents Worldwide:	Number of potential contributing factors: Yes - some potential factors (up to 3)			Reusable: No	Similar events (past 6 months): 1 incidents
Products supplied the last 12 months:	Specific factors identified: Compatibility of device - patient characteristics, Use of device - operator/user error	ESTIMATED LEVEL OF PRIORITY: Routine	FINAL LEVEL OF PRIORITY: Routine		3 or more events - batch/model: No
Products supplied last 24 months:	Number of potential sensitivities: No	EXCEPTION TO PRIORITY LEVEL: Low incidence			3 or more events - health district: No
Products supplied last 36 months:	Specific sensitivities identified:				3 or more events - organisation: No
Products supplied Worldwide:	Consultations during risk assessment: I undertook an internet search (e.g., Google)	Final Risk Assessment: No			

Additional Risk Analysis


Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section					
Update Device Details?:	Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
Yes	Date: 24/02/2020	Severity: 5 - An illness/injury was resolved or prevented with treatment by a health professional	Incidents in the last 12 months:	Manufacturer analysis: Yes		
Copy Data From: Initial						
	Assessor: §22	Manufacturer documentation: Unknown - updated information from the manufacturer is required	Incidents in last 24 months:	Manufacturer action: No	ESTIMATED LEVEL OF INVESTIGATION: Level 1 Investigation (to complete screening)	FINAL LEVEL OF INVESTIGATION: Screening only
	Injured Party: Patient	Device Recalls: 0. No recalls for similar incidents in Australia	Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL: Closed with NFA. Device not returned.	
	Found Prior To Use: No	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors: Yes - some potential factors (up to 3)		
	Reusable: No	Similar events (past 6 months): 0 incidents	Products supplied the last 12 months:	Specific factors identified: Compatibility of device - patient characteristics, Use of device - operator/user error	ESTIMATED LEVEL OF PRIORITY: Routine	FINAL LEVEL OF PRIORITY: Routine
		3 or more events - batch/model: No	Products supplied last 24 months:	Number of potential sensitivities: No	EXCEPTION TO PRIORITY LEVEL:	
		3 or more events - health district: No	Products supplied last 36 months:	Specific sensitivities identified:		
		3 or more events - organisation: No	Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment: Yes	

Sponsor/Manufacturer Information Section

<b>Search Sponsors:</b> 837	<b>Name:</b> Medtronic Australasia Pty Ltd	<b>Client #:</b> 837
<b>Attention To:</b> §22	<b>Address 1:</b> PO Box 945	<b>Address 2:</b> 
<b>State:</b> NSW	<b>Postcode:</b> 1670	<b>Town/Suburb:</b> NORTH RYDE BC
<b>Email:</b> §22@medtronic.com	<b>Phone:</b> §22	<b>Fax:</b> 

Investigation Information Section - Submitted by Sponsor/Manufacturer

<b>Device Analysis Results:</b> The product sample or additional supporting materials from the account was not available for this incident. Records from each manufacturing lot are thoroughly reviewed to ensure that products are released meeting all Medtronic quality release specifications at the time of manufacture. There were no assembly or component issues identified; therefore, a Device History Record review will not be performed. Without the physical product, a definitive root cause could not be identified. Post Market Vigilance will continue to monitor this condition for future potential action. Should new information become available, the file will be re-opened and the investigation will be amended as appropriate.	<b>Corrective/Preventative Actions:</b> None at this time
<b>Details of Similar Events:</b> 0 Similar events (0%) from 2617 sales in Australia (past 3 years) . 12 Similar events (0.037%) from 32248 sales worldwide (past 3 years).	<b>Additional Details (use for tables):</b> 
<b>CAPA# Reference:</b> 	
<b>Risk Assessment</b>	
<b>Frequency:</b> 	<b>Severity:</b> 
<b>Rating:</b> 	<b>Type Cause and Outcome:</b> 
<b>Expected Rate:</b> 	<b>Number of Similar Events:</b> 
<b>Actual Rate:</b> 	
<b>Countries Similar Events Also Occurred:</b> United States, United Kingdom, Netherlands, France, Australia, Sweden.	
<b>Completed Actions:</b> 	<b>Planned Actions and Proposed Timelines:</b> 
<b>Additional Comments:</b> DIR closed 4/3/20	

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes
<input type="checkbox"/>		TGA Regulatory Action	Final Report Due		26/11/2019	05/03/2020	24/02/2020		

List of Problem Observed Codes - Click [N] to begin entering information.

<b>Problem Observed Details</b>	
<b>Problem Observed (Level 1)</b>	<b>Problem Observed (Level 2)</b>
<b>Problem Observed (Level 3)</b>	<b>If 'Other' Selected</b>

Appropriate Term/Code Not Available			anastomosis dehiscence	
Appropriate Term/Code Not Available			pelvic abscess	
Appropriate Term/Code Not Available			tissue loss	

Clinical signs symptoms and conditions

Details			
Level 1	Level 2	Level 3	
Procedural Complications	Wound Dehiscence		
Infections	Abscess		

Health Impact

Details			
Level 1	Level 2	Level 3	
Surgical Intervention	Additional Surgery		

Investigation Findings

Finding Details			
Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
No Findings Available			

Investigation Conclusion

Conclusion Details			
Investigation Conclusion (L1)	Investigation Conclusion (L2)	If Additional Conclusion Detail Requested	
Cause Not Established			

Investigation Outcomes

Outcome Details			
Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested	
Reviewed, for Trending Purposes Only			

Investigation Summary

Investigation Type: <input type="text"/>	Latest Investigation (DII) where this DIR is the Primary DIR: <input type="text"/>	Latest Investigation (DII) where this DIR is a Related DIR: <input type="text"/>	Investigator: <input type="text"/>	Extension Number: <input type="text"/>
Investigator's Notes: <input type="text"/>		Summary Findings: <p>No further investigation will occur at this time, however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.</p>		Recall Number: <input type="text"/>

**Note:** Letter generation buttons disabled if report not ready for website publication or risk analysis not completed.

This section is used to match information provided via UDIR forms to ARTG information. You can select a Brand/Name from information provided in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered): <input type="text"/>	Brand Name: <input type="text"/>	Manufacturer Name: <input type="text"/>	Device ARTG #: <input type="text"/>
---	-------------------------------------	--	--

Other Devices					
Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Click **New** to begin entering information.

Rec No	
1	

Samples Record - Click **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details			Additional Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Device Description:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Version Number:	
					Who sent the device to the TGA?:			Why does the TGA have the sample?:	

**Additional Patients**

Click **[N]** to begin entering information.

Patient Details			
Sex:	Weight:	Age:	
Patient Focused Corrective Action Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

**Additional Device Information**

Where did you get this device from?:  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

Similar Events

Similar events - how many times?:

Date of Recent Report:

Event Reported To:

Reporter Reference Number:

Device Access - Alternate Device Contact Information Provided

Title:

First Name:

Fax:

Email:

Last Name:

Phone:

Incident Location Details

Occurred in Australia:

Organisation:

Town/Suburb:

State:

Address Line 1:

Address Line 2:

Postcode:

Flow Details : DIR-REQ - Device Incident Request : 220343

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
220343	DIR-REQ		Closed	522	OPR Administration User	04/03/2020	Normal	0

Signature Details

Role	IRIS Investigator	
User	522	
Signed At	13/08/2020 11:13:35	
Comment		



## Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 38 - ID : 493698

Released by **522** on 14/04/2020 09:01:24

Report #:  Records Management #:  Reporter's Reference #:  Report Type:

ARTG: 293808

[Document Container URL](#)

### Report Information Section

Report Status:  Sponsor's Reported Category:  Date of Adverse Event:  Date of Initial Report:

Date of Final Report:  Date of Initial TGA Act on:  Reviewed by Team:  Date Response Received:

Date Completed:  Operator at Time of Event:  If 'Other' Operator Selected:  Reporter consents to contact by sponsor:

Source of Report:  If 'Other' Source Selected:  Type of Initial Action:

### Event Description for Website Publication:

### Clinical Event Information:

Number of Incidents in Report:  Contact:  Alternative Person Title:  Alternative Person First Name:

Alternative Person Surname:  Alternative Person Phone:  Alternative Person Fax:  Alternative Person Email:

### Recorded Problems Observed

#### Recorded Problems Observed:

### Clinical Signs, Symptoms and Conditions

#### Recorded Clinical Signs, Symptoms and Conditions:

### Health Impact

#### Recorded Health Impacts:

### Patient Information

Sex:  Weight:  Age:

Patient Focused Corrective Action Taken:  Patient History:

Patient Outcome/Consequences:  Additional Event Description:

Describe any test (Lab, xray, etc.):  Injured - Extent of Injury:  Other medical devices currently using/implanted:  Medical Problem Device Used For:

### Additional Patients Added:



ADD OTHER PATIENTS ABOVE.

0

Submitting Reporter Section

Search Reporter By Surname:

Reporter Title:

Position:

Address 1:

Country:

Mobile:

Reporter #:

First Name:

Surname:

Company/Inst tution:

Town/Suburb:

State:

Phone:

Fax:

Last External Submission By:

Initial Reporter Section

As Above?:

If No, fill out the following:

Initial Reporter Confidential:

Search Reporter By Surname:

Initial Reporter #:

Preferred Contact Method:

Title:

First Name:

Surname:

Position:

Company/Inst tution:

Address 1:

Address 2:

Town/Suburb:

State:

Postcode:

Country:

Phone:

Fax:

Mobile:

Email:

Allow the device company to contact you about the incident:

Device Information Section

Product Exempt (Note: If not exempt, enter ARTG No):

Search Device ARTG:

Device ARTG #:

Therapeutic Licence Type:

Product Licence Category:

Device Class:

GMDN / UMDN Code:

GMDN / UMDN Text:

Brand Name:

Initial Device Description:

Usage of Device:

Software Version:

Model #:

Serial #:

Batch #:

Lot #:

Purchase Date:

Expiry Date:

Date of Implant:

Date of Explant:

Date of Initial Procedure:

Place of Implantation:

Reported Device Location:

Access Contact Title:

Access Contact First Name:

Access Contact Surname:

Access Contact Phone:

Access Contact Fax:

Access Contact Email: <input type="text"/>	Licence Status: <input type="text" value="A"/>	Status Effective Date: <input type="text" value="13/09/2017"/>	Additional Devices Added: <input type="text" value="0"/>
---	---	---	---

Manufacturer Information Section

Manufacturer Name: <input type="text" value="Covidien Inc"/>	Manufacturer Client Id: <input type="text" value="54968"/>	Address 1: <input type="text"/>
Address 2: <input type="text"/>	State/Province: <input type="text"/>	Country: <input type="text"/>
Postcode: <input type="text"/>	Phone: <input type="text"/>	Email: <input type="text"/>
Manufacturer Informed: <input type="text" value="Yes"/>	Date Aware of Adverse Event: <input type="text" value="25/03/2020"/>	Contact First Name: <input type="text"/>
Contact Surname: <input type="text"/>	Contact Title: <input type="text"/>	

Supplier Information Section

Supplier Name: <input type="text"/>	Address 1: <input type="text"/>	Address 2: <input type="text"/>
Town/Suburb: <input type="text"/>	Country: <input type="text"/>	Postcode: <input type="text"/>
Phone: <input type="text"/>	Email: <input type="text"/>	Website: <input type="text"/>
Supplier Informed: <input type="text"/>	Contact Title: <input type="text"/>	Contact First Name: <input type="text"/>
Contact Surname: <input type="text"/>	Contact Fax: <input type="text"/>	Contact Email: <input type="text"/>

Report Status					
For website publication: <input type="text" value="Yes"/>	Ready for Publication: <input type="text" value="Yes"/>	Investigated: <input type="text" value="No"/>	Investigation Reason: <input type="text" value="Event determined to be most likely due to user error"/>	Team Assignment: <input type="text" value="Team B (IIB, IS &amp; IM)"/>	Team Priority: <input type="text" value="Not Investigated"/>

Team Review		
Reviewed by Team: <input type="text"/>	Reason Sent To Meeting: <input type="text"/>	Outcome from team meeting: <input type="text"/>
Notes for Team meeting: <input type="text"/>		
Outcomes from Team Meeting: <input type="text"/>		

DPRC Review		
Reviewed by DPRC: <input type="text"/>	DPRC Reason Sent To Meeting: <input type="text"/>	Outcome from DPRC Meeting: <input type="text"/>
Meeting Notes: <input type="text"/>		

Initial Risk Analysis					
Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	Assessor:	Manufacturer documentation:

27/03/2020	1 - No harm has been reported to occur			Chloe Chen	Unknown - updated informat on from the manufacturer is required
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:	Injured Party:	Device Recalls:
		Level 1 Investigat on (to complete screening)	Screening only	Not Appl cable	0. No recalls for similar inc dents in Australia
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:		Found Prior To Use:	Is AE covered by current recall:
				Yes	
Incidents Worldw de:	Number of potential contributing factors:			Reusable:	Similar events (past 6 months):
	No			No	0 inc dents
Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:	FINAL LEVEL OF PRIORITY:	3 or more events - batch/model:	
		Routine	Routine		
Products supplied last 24 months:	Number of potential sens tivities:	EXCEPTION TO PRIORITY LEVEL:		3 or more events - health district:	
	No				
Products supplied last 36 months:	Specific sens tivities identified:			3 or more events - organisation:	
Products supplied Worldw de:	Consultations during risk assessment:	Final Risk Assessment:			
		No			

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section				
Update Device Details?:	<b>Background Information</b>	<b>Risk Assessment - Section A</b>	<b>Risk Assessment - Section B</b>	<b>Risk Assessment - Section C</b>	<b>Risk Assessment - Section D</b>
Yes	Date:	Severity:	Inc dents in the last 12 months:	Manufacturer analysis:	
Copy Data From:	05/08/2020	<b>5 - An illness/injury was resolved or prevented with treatment by a health professional</b>		Yes	
In tial					
	Assessor:	Manufacturer documentation:	Inc dents in last 24 months:	Manufacturer act on:	ESTIMATED LEVEL OF INVESTIGATION:
	S22	<b>Unknown - updated information from the manufacturer is required</b>		Not required (dev ce was performing within specified parameters)	Level 1 Investigat on (to complete screening)
	Injured Party:	Device Recalls:	Inc dents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:
	<b>Patient</b>	<b>0. No recalls for similar incidents in Australia</b>			Closed. no recall or PMRs. Low rates qlik and sponsor analysis.
	Found Pr or To Use:	Is AE covered by current recall:	Inc dents Worldw de:	Number of potential contributing factors:	Likely to be user error.
	<b>Yes</b>			Yes - some potential factors (up to 3)	
	Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:
	<b>No</b>	<b>0 incidents</b>		Use of device - operator/user error	Routine
		3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitiv ties:	EXCEPTION TO PRIORITY LEVEL:
		No		Yes - some potential sensitiv ties (up to 3)	Closed. no recall or PMRs. Low rates qlik and sponsor analysis.
		3 or more events - health district:	Products supplied last 36 months:	Specific sensitiv ties identified:	Likely to be user error.
		No		Dev ce used in high risk populations, Dev ce used in high acu ty clinical environments	
		3 or more events - organisation:	Products supplied Worldwide:	Consultat ons during risk assessment:	Final Risk Assessment:
		No		I did none of the above incidents	<b>Yes</b>


Sponsor/Manufacturer Information Section

Search Sponsors:	Name:	Client #:
------------------	-------	-----------

837	Medtron c Australasia Pty Ltd	837
Attention To:	Address 1:	Town/Suburb:
§22	PO Box 945	NORTH RYDE BC
State:	Postcode:	Phone:
NSW	1670	§22
Email:		Fax:
§22 @medtron c.com		

Investigat on Information Sect on - Submitted by Sponsor/Manufacturer

<p><b>Device Analysis Results:</b></p> <p>Medtronic conducted an investigat on based upon all received informat on.</p> <p>The following complaint(s) were investigated: It was reported that: Significant surgical intervent on was required due to the product failure Significant tissue loss occurred as a result of the product failure The anvil tilted prematurely / pr or to firing. - This complaint was confirmed based on the actual device. The most likely root cause was traced to the user. The anvil did not easily attach to the stapler - This complaint was confirmed based on the actual device. The most likely root cause was traced to the user.</p> <p>A medical assessment was performed. The medical assessment determined there was no association to the complaint. A risk assessment review was performed. The risk assessment review determined that there was a relat on to the complaint.</p> <p>Product analysis occurred. The primary analysis finding was: - Damage consistent with an obstruction was noted on the knife blade and/or staple guide.</p>	<p><b>Corrective/Preventative Actions:</b></p> <p>Further action was not required because the event had foreseen risk and is included in a data mon toring plan. Should new informat on become available the file will be re-opened and the investigat on summary will be amended as appropriate.</p>
---	---

<p><b>Details of Similar Events:</b></p> <p>Similar inc dents based off FDA Annex A code A040102 (LOSS OF OR FAILURE TO BOND): Similar events in Australia in the past 3 years: 2 Similar events worldw de (Australia excluded) in the past 3 years: 22 Sales in Australia in the past 3 years: 4376 This correlates to a complaint rate of 0.046% Sales worldwide (Australia excluded) in the past 3 years: 61876 This correlates to a complaint rate of 0.036%</p> <p>Similar inc dents based off FDA Annex A code A051204 (DEVICE SLIPPED): Similar events in Australia in the past 3 years: 1 Similar events worldw de (Australia excluded) in the past 3 years: 1 Sales in Australia in the past 3 years: 4376 This correlates to a complaint rate of 0.023% Sales worldwide (Australia excluded) in the past 3 years: 61876 This correlates to a complaint rate of 0.002%</p>	<p><b>Additional Details (use for tables):</b></p> 
---	--

CAPA# Reference:

<b>Risk Assessment</b>		<b>Type Cause and Outcome:</b>	<b>Number of Similar Events:</b>
<b>Frequency:</b>	<b>Sever ty:</b>		
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<b>Rating:</b>			
<input type="text"/>			
<b>Expected Rate:</b>	<b>Actual Rate:</b>		
<input type="text"/>	<input type="text"/>		

**Countries Similar Events Also Occurred:**

<b>Completed Actions:</b>	<b>Planned Act ons and Proposed Timelines:</b>
<input type="text"/>	<input type="text"/>

**Addit onal Comments:**

Details of Reasons

Reason for Level 1 Investigat on

Focus of Level 2 Investigation

Details of Focus

Essential Principles

If 'Other' Selected

Sources of Evidence for Level 2

Details of Source

Sources of Ev dence

If 'Others' please specify here

Expected Sourcing Date

Date of Ev dence Received

Evidence

Investigat on Quest ons (Level 1 and Level 2):

Potential Risks

Delays in response by product manufacturers:

Delays in response by incident reporters:

Delays in analysis w thin the TGA:

Delays in reporting by other sources (e.g. clin cal registries):

Other Risks (wh ch need to be specified):

Next Steps for Level 1 & Level 2 Investigat ons

Next Steps for Level 1 Investigat on:

Next Steps for Level 2 Investigat on:

Cl ck [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes
<input type="checkbox"/>		TGA Regulatory Action	Final Report Due		30/03/2020	08/07/2020	31/07/2020		

List of Problem Observed Codes - Cl ck [N] to begin entering informat on.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Compatibility Problem	Component or Accessory Incompatibility		
Mechan cal Problem	Mechan cal Jam		

Clinical signs symptoms and conditions

Details

Level 1	Level 2	Level 3	
Others	Appropriate Term / Code Not Available		

Health Impact

Details			
Level 1	Level 2	Level 3	
Unexpected Medical Intervention	Additional Device Required		

Investigation Findings

Finding Details			
Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
Mechanical Problem Identified			likely due to User Error

Investigation Conclusion

Conclusion Details		
Investigation Conclusion (L1)	Investigation Conclusion (L2)	If Additional Conclusion Detail Requested
Cause Traced to User		

Investigation Outcomes

Outcome Details		
Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested
Reviewed, for Trending Purposes Only		

Investigation Summary

Investigation Type:	Latest Investigation (DII) where this DIR is the Primary DIR:	Latest Investigation (DII) where this DIR is a Related DIR:	Investigator:	Extens on Number:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Investigator's Notes:	Summary Findings:		Recall Number:	
<input type="text"/>	No further investigation will occur at this time; however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.		<input type="text"/>	

**Note:** Letter generation buttons disabled if report not ready for website publication or risk analysis not completed.

Device Lookup

This section is used to match information provided via UDIR forms to ARTG information. You can select a Brand/Name from information provided in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered):	Brand Name:	Manufacturer Name:	Device ARTG #:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Other Devices					
Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:

Model Number:	Batch #:	Lot #:	Expiry Date:	
---------------	----------	--------	--------------	--

Related DIR Information - Click **New** to begin entering information.

Rec No	
1	

Samples Record - Click **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details			Additional Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Device Description:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Version Number:	
					Who sent the device to the TGA?:			Why does the TGA have the sample?:	

**Additional Patients**

Click **[N]** to begin entering information.

<b>Patient Details</b>			
Sex:	Weight:	Age:	
Patient Focused Corrective Act on Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

**Additional Device Information**

Where did you get this device from?:  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

**Similar Events**

Similar events - how many times?:  Date of Recent Report:  Event Reported To:  Reporter Reference Number:

**Device Access - Alternate Device Contact Information Provided**

Title:  First Name:  Last Name:  Phone:

Fax:  Email:

Two empty input fields.

Incident Location Details

Occurred in Australia:

Empty input field.

Organisation:

Empty input field.

Town/Suburb:

Empty input field.

State:

Empty input field.

Address Line 1:

Empty input field.

Address Line 2:

Empty input field.

Postcode:

Empty input field.

Flow Details DIR-REQ - Device Incident Request 259093

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
259093	DIR-REQ		Closed	S22	OPR Administration User	05/08/2020	Normal	0

Signature Details

Role	IRIS Investigator	
User	S22	
Signed At	05/08/2020 13:59:18	
Comment		





# Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 44 - ID : 528946

Released by 522 on 18/03/2021 17:26:54

Report #: <input type="text" value="72184"/>	Records Management #: <input type="text"/>	Reporter's Reference #: <input type="text" value="704585197"/>	Report Type: <input type="text" value="Final"/>
---	---	---	--

ARTG: 293808	<a href="#">Document Container URL</a>		
--------------	--	--	--

### Report Information on Section

Report Status: <input type="text" value="Closed"/>	Sponsor's Reported Category: <input type="text" value="Other"/>	Date of Adverse Event: <input type="text"/>	Date of Initial Report: <input type="text" value="03/09/2021"/>
Date of Final Report: <input type="text" value="21/10/2021"/>	Date of Initial TGA Action: <input type="text" value="03/09/2021"/>	Reviewed by Team: <input type="text"/>	Date Response Received: <input type="text"/>
Date Completed: <input type="text" value="25/10/2021"/>	Operator at Time of Event: <input type="text"/>	If 'Other' Operator Selected: <input type="text"/>	Reporter consents to contact by sponsor: <input type="text" value="N/A"/>
Source of Report: <input type="text" value="Sponsor"/>	If 'Other' Source Selected: <input type="text"/>	Type of Initial Act on: <input type="text" value="Trend data only"/>	

Event Description for Website Publication:

Clinical Event Information:

Number of Incidents in Report: <input type="text" value="1"/>	Contact: <input type="text"/>	Alternative Person Title: <input type="text"/>	Alternative Person First Name: <input type="text"/>
Alternative Person Surname: <input type="text"/>	Alternative Person Phone: <input type="text"/>	Alternative Person Fax: <input type="text"/>	Alternative Person Email: <input type="text"/>

### Recorded Problems Observed

Recorded Problems Observed:

### Clinical Signs, Symptoms and Conditions

Recorded Clinical Signs, Symptoms and Conditions:

### Health Impact

Recorded Health Impacts:

### Patient Information

Sex: <input type="text"/>	Weight: <input type="text"/>	Age: <input type="text"/>	
Patient Focused Corrective Action Taken: <input type="text"/>		Patient History: <input type="text"/>	
Patient Outcome/Consequences: <input type="text"/>		Additional Event Description: <input type="text"/>	
Describe any test (Lab, xray, etc.): <input type="text"/>	Injured - Extent of Injury: <input type="text"/>	Other medical devices currently using/implanted: <input type="text"/>	Medical Problem Device Used For: <input type="text"/>

Additional Patients Added:

0

Submitting Reporter Section

<b>Search Reporter By Surname:</b> §22	<b>Reporter #:</b> 	<b>Preferred Contact Method:</b> 
<b>Reporter Title:</b> §22	<b>First Name:</b> §22	<b>Surname:</b> §22
<b>Position:</b> §22	<b>Company/Instituti on:</b> Medtronic Australasia Pty Ltd	<b>State:</b> NSW
<b>Address 1:</b> 2 Alma Road	<b>Address 2:</b> 	<b>Town/Suburb:</b> Macquarie Park
<b>Country:</b> Australia	<b>Postcode:</b> 2113	<b>Phone:</b> §22
<b>Mobile:</b> 	<b>Email:</b> §22 @medtron c.com	<b>Last External Submissi on By:</b> 113291_837 - 21/10/2021 15:42

Initial Reporter Section

<b>As Above?:</b> No	<i>If No, fill out the following:</i>	<b>Initial Reporter Confidential:</b> Yes
<b>Search Reporter By Surname:</b> 	<b>Initial Reporter #:</b> 	<b>Preferred Contact Method:</b> 
<b>Title:</b> 	<b>First Name:</b> 	<b>Surname:</b> 
<b>Position:</b> 	<b>Company/Instituti on:</b> KNOX PRIVATE HOSPITAL	<b>State:</b> VIC
<b>Address 1:</b> 262 Mountain Highway	<b>Address 2:</b> 	<b>Town/Suburb:</b> Wantirna
<b>Postcode:</b> 3152	<b>Country:</b> 	<b>Phone:</b> §22
<b>Mobile:</b> 	<b>Email:</b> ap.callcentre@healthscope.com	<b>Allow the device company to contact you about the incident:</b> <input type="checkbox"/>

Device Information Section

<b>Product Exempt (Note: If not exempt, enter ARTG No):</b> No	<b>Search Device ARTG:</b> 293808	<b>Device ARTG #:</b> 293808	<b>Therapeutic Licence Type:</b> Medical Device
<b>Product Licence Category:</b> Included	<b>Device Class:</b> Class IIb	<b>GMDN / UMDN Code:</b> 59875	<b>GMDN / UMDN Text:</b> Intraluminal circular stapler, single-use
<b>Brand Name:</b> EEA - Intraluminal circular stapler, single-use	<b>Initial Device Description:</b> EEA - Intraluminal circular stapler, single-use	<b>Usage of Device:</b> Single Use	<b>Software Versions:</b> 
<b>Model #:</b> TRIEEA28MT	<b>Serial #:</b> Unknown	<b>Batch #:</b> 	<b>Lot #:</b> 
<b>Purchase Date:</b> 	<b>Expiry Date:</b> 	<b>Date of Implant:</b> 	<b>Date of Explant:</b> 
<b>Date of Initial Procedure:</b> 	<b>Place of Implantation:</b> 	<b>Reported Device Location:</b> 	<b>Access Contact Title:</b> 
<b>Access Contact First Name:</b> 	<b>Access Contact Surname:</b> 	<b>Access Contact Phone:</b> 	<b>Access Contact Fax:</b> 
<b>Access Contact Email:</b> 	<b>Licence Status:</b> 	<b>Status Effective Date:</b> 	<b>Additional Devices Added:</b> 

<input type="text"/>	<input type="text" value="A"/>	<input type="text" value="13/09/2017"/>	<input type="text" value="0"/>
Manufacturer Information Section			
<b>Manufacturer Name:</b> <input type="text" value="Covidien llc"/>	<b>Manufacturer Client Id:</b> <input type="text" value="54968"/>	<b>Address 1:</b> <input type="text"/>	
<b>Address 2:</b> <input type="text"/>	<b>Town/Suburb:</b> <input type="text"/>	<b>State/Province:</b> <input type="text"/>	<b>Country:</b> <input type="text"/>
<b>Postcode:</b> <input type="text"/>	<b>Phone:</b> <input type="text"/>	<b>Fax:</b> <input type="text"/>	<b>Email:</b> <input type="text"/>
<b>Manufacturer Informed:</b> <input type="text" value="Yes"/>	<b>Date Aware of Adverse Event:</b> <input type="text" value="31/08/2021"/>	<b>Contact Title:</b> <input type="text"/>	<b>Contact First Name:</b> <input type="text"/>
<b>Contact Surname:</b> <input type="text"/>			

Supplier Information Section

<b>Supplier Name:</b> <input type="text"/>	<b>Address 1:</b> <input type="text"/>	<b>Address 2:</b> <input type="text"/>
<b>Town/Suburb:</b> <input type="text"/>	<b>State:</b> <input type="text"/>	<b>Postcode:</b> <input type="text"/>
<b>Phone:</b> <input type="text"/>	<b>Country:</b> <input type="text"/>	<b>Website:</b> <input type="text"/>
<b>Supplier Informed:</b> <input type="text"/>	<b>Date of Supplier Contact:</b> <input type="text"/>	<b>Contact First Name:</b> <input type="text"/>
<b>Contact Surname:</b> <input type="text"/>	<b>Contact Phone:</b> <input type="text"/>	<b>Contact Email:</b> <input type="text"/>
	<b>Contact Title:</b> <input type="text"/>	
	<b>Contact Fax:</b> <input type="text"/>	

Report Status

<b>For website publication:</b> <input type="text" value="Yes"/>	<b>Ready for Publication:</b> <input type="text" value="Yes"/>	<b>Investigated:</b> <input type="text" value="No"/>	<b>Investigation Reason:</b> <input type="text" value="Rate considered low at this stage"/>	<b>Team Assignment:</b> <input type="text" value="Team B (IIB, IS &amp; IM)"/>	<b>Team Priority:</b> <input type="text" value="Not Investigated"/>
---	---	---	--	---	--

Team Review

<b>Reviewed by Team:</b> <input type="text"/>	<b>Reason Sent To Meeting:</b> <input type="text"/>	<b>Outcome from team meeting:</b> <input type="text"/>	<input type="text"/>
<b>Notes for Team meeting:</b> <input type="text"/>			
<b>Outcomes from Team Meeting:</b> <input type="text"/>			

Initial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
<b>Date:</b> <input type="text" value="03/09/2021"/>	<b>Severity:</b> <input type="text" value="5 - An illness/injury was resolved or prevented with treatment by a health professional"/>	<b>Incidents in the last 12 months:</b> <input type="text"/>	<b>Manufacturer analysis:</b> <input type="text"/>	<b>Assessor:</b> <input type="text" value="S22"/>	<b>Manufacturer documentation:</b> <input type="text"/>
<b>Incidents in last 24 months:</b> <input type="text"/>	<b>Manufacturer action:</b> <input type="text"/>	<b>ESTIMATED LEVEL OF INVESTIGATION:</b> <input type="text" value="Screening only"/>	<b>FINAL LEVEL OF INVESTIGATION:</b> <input type="text" value="Screening only"/>	<b>Injured Party:</b> <input type="text" value="Patient"/>	<b>Device Recalls:</b> <input type="text"/>
<b>Incidents in last 36 months:</b> <input type="text"/>	<b>IVD status:</b> <input type="text"/>	<b>EXCEPTION TO INVESTIGATION LEVEL:</b> <input type="text"/>		<b>Found Prior To Use:</b> <input type="text" value="No"/>	<b>Is AE covered by current recall:</b> <input type="text"/>
<b>Incidents Worldwide:</b>	<b>Number of potential contributing factors:</b>			<b>Reusable:</b>	<b>Similar events (past 6 months):</b>

<input type="text"/>	No	<input type="text"/>	No	<input type="text"/>
Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:	FINAL LEVEL OF PRIORITY:	0 incidents
<input type="text"/>	<input type="text"/>	Routine	Routine	3 or more events - batch/model:
Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:		3 or more events - health district:
<input type="text"/>	No	<input type="text"/>		3 or more events - organisation:
Products supplied last 36 months:	Specific sensitivities identified:			<input type="text"/>
<input type="text"/>	<input type="text"/>			
Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:		
<input type="text"/>	I did none of the above incidents	No		

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section				
Update Device Details?:	Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
Yes	Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	
Copy Data From:	25/10/2021	5 - An illness/injury was resolved or prevented with treatment by a health professional		No	
In trial					
	Assessor:	Manufacturer documentation:	Incidents in last 24 months:	Manufacturer act on:	ESTIMATED LEVEL OF INVESTIGATION:
	S22	Yes - manufacturer HAS documented issues or complications of this type		Unable to take action (device was not returned for analysis)	Screening only
	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:			Screening only
	Injured Party:	Device Recalls:	Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:
	Patient	0. No recalls for similar incidents in Australia			
	Found Prior To Use:	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors:	
	No	No		No	
	Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:
	No	0 incidents			Routine
		3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitivities:	FINAL LEVEL OF PRIORITY:
				No	Routine
		3 or more events - health district:	Products supplied last 36 months:	Specific sensitivities identified:	EXCEPTION TO PRIORITY LEVEL:
		3 or more events - organisation:	Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:
					Yes

Sponsor/Manufacturer Information Section

Search Sponsors:	Name:	Client #:
837	Medtronic Australasia Pty Ltd	837
Attention To:	Address 1:	Town/Suburb:
S22	PO Box 945	NORTH RYDE BC
State:	Postcode:	Phone:
NSW	1670	S22
Email:		Fax:
S22@medtronic.com		

Investigat on Information Sect on - Submitted by Sponsor/Manufacturer

Device Analysis Results:

Medtronic conducted an investigat on based upon all received informat on.

The following complaint(s) were investigated:

It was reported that:

An unexpected content leak occurred along the staple line.

- There was not enough informat on to make any determination based on insuff cient evidence. The most likely root cause was not established.

The deployed staples did not hold the tissue and the staple line opened up

- There was not enough informat on to make any determination based on insuff cient evidence. The most likely root cause was not established.

A device history review was performed and the lot number was not prov ded. However, records from each manufacturing lot are thoroughly reviewed to ensure that products are released meeting all Medtron c quality release specificat ons at the time of manufacture. . The device history review determined that there was insuff cient informat on to conduct the specific Investigat on Action.

A medical assessment was performed. The medical assessment results were inconclusive.

Further action was not required because there was insufficient information.

Should new informat on become available the file will be re-opened and the investigation summary will be amended as appropriate.

Details of Similar Events:

Australia:

0 similar events from 3993 units (0%) sold in CY2019;

0 similar events from 2977 units (0%) sold in CY2020;

1 similar events from 2744 units (0.03644%) sold in CY2021 (Includes this event).

Worldw de:

8 similar events from 46210 un ts (0.01731%) sold in CY2019;

8 similar events from 41670 un ts (0.0192%) sold in CY2020;

21 similar events from 46093 units (0.04556%) sold in CY2021.

Similar inc dent rates were determined using FDD annex A code (ADVERSE EVENT WITHOUT IDENTIFIED DEVICE OR USE PROBLEM - A24) for the platform family type of product reported per subject event.

Addit onal Details (use for tables):



CAPA# Reference:

Risk Assessment

Frequency:

Severity:

Rating:

Type Cause and Outcome:

Number of Similar Events:

Expected Rate:

Actual Rate:

Countries Similar Events Also Occurred:

Albania Argentina Armenia Australia Austria Bahrain Belgium Brazil Bulgaria Canada Canary Islands Chile China Colombia Costa Rica Croatia Cyprus Czech Republ c Denmark Finland France Germany Greece Guadeloupe Guatemala Hong Kong Hungary Iceland India Indonesia Ireland Israel Italy Korea, Republ c Of Kuwa t Latvia Lebanon Luxembourg Macao Malaysia Montenegro Netherlands New Zealand North Macedonia Norway Pakistan Panama Paraguay Peru Poland Portugal Puerto R co Qatar Reunion RomaniaRussian Federation San Marino Saudi Arabia Serbia Singapore Slovakia Slovenia South Afr ca Spain Sweden Sw tzerland Thailand Turkey UkraineUn ted Arab Emirates Un ted Kingdom Un ted States Viet Nam

Completed Actions:

Planned Actions and Proposed Timelines:

Addit onal Comments:

DIR closed on 25 OCT 2021

Should new informat on become available the file will be re-opened and the investigation summary will be amended as appropriate.

Reason for Level 1 Investigation

Details of Reasons

Reason for Level 1 Investigat on

Focus of Level 2 Investigation

Details of Focus

Essential Principles

If 'Other' Selected

Sources of Evidence for Level 2

Details of Source			
Sources of Evidence	If 'Others' please specify here	Expected Sourcing Date	Date of Evidence Received

Evidence

Investigation Questions (Level 1 and Level 2):

Potential Risks

Delays in response by product manufacturers:

Delays in response by incident reporters:

Delays in analysis within the TGA:

Delays in reporting by other sources (e.g. clinical registries):

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:

Next Steps for Level 2 Investigation:

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Insufficient Information			

Clinical signs symptoms and conditions

Details			
Level 1	Level 2	Level 3	
Gastrointestinal System			

Health Impact

Details			
Level 1	Level 2	Level 3	
Surgical Intervention	Additional Surgery		

Investigation Findings

Finding Details	



Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Vers on Number:	<b>Document 4</b>
				Who sent the device to the TGA?:			Why does the TGA have the sample?:	

**Additional Patients**

Click **[N]** to begin entering information.

<b>Patient Details</b>			
Sex:	Weight:	Age:	
Patient Focused Corrective Act on Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

**Additional Device Information**

Where did you get this device from?:  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

**Similar Events**

Similar events - how many times?: <input type="text"/>	Date of Recent Report: <input type="text"/>	Event Reported To: <input type="text"/>	Reporter Reference Number: <input type="text"/>
--	---	---	---

**Device Access - Alternate Device Contact Information Provided**

Title: <input type="text"/>	First Name: <input type="text"/>	Last Name: <input type="text"/>	Phone: <input type="text"/>
Fax: <input type="text"/>	Email: <input type="text"/>		

**Incident Location Details**

Occurred in Australia: <input type="text"/>	Organisation: <input type="text"/>	Address Line 1: <input type="text"/>	Address Line 2: <input type="text"/>
Town/Suburb: <input type="text"/>	State: <input type="text"/>	Postcode: <input type="text"/>	



ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
318468	DIR-REQ		Closed	<b>s22</b>	OPR Administration User	25/10/2021 ✓	Normal	0

Signature Details

Role	IRIS Investigator	
User	<b>s22</b>	
Signed At	25/10/2021 19:04:18	
Comment		



# Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 45 - ID : 538597

Released by Theta Technologies on 24/11/2021 14:57:03

Report #:	Records Management #:	Reporter's Reference #:	Report Type:
<input type="text" value="76781"/>	<input type="text"/>	<input type="text" value="0704894000"/>	<input type="text" value="Final"/>

ARTG: 293808	<a href="#">Document Container URL</a>
--------------	--

### Report Information Section

Report Status:	Sponsor's Reported Category:	Date of Adverse Event:	Date of Initial Report:
<input type="text" value="Closed"/>	<input type="text"/>	<input type="text" value="S22"/>	<input type="text" value="08/03/2022"/>
Date of Final Report:	Date of Initial TGA Action:	Reviewed by Team:	Date Response Received:
<input type="text" value="07/07/2022"/>	<input type="text" value="08/03/2022"/>	<input type="text" value="19/07/2022"/>	<input type="text"/>
Date Completed:	Operator at Time of Event:	If 'Other' Operator Selected:	Reporter consents to contact by sponsor:
<input type="text" value="19/07/2022"/>	<input type="text" value="Healthcare Professional"/>	<input type="text"/>	<input type="text" value="N/A"/>
Source of Report:	If 'Other' Source Selected:	Type of Initial Action:	
<input type="text" value="Sponsor"/>	<input type="text"/>	<input type="text" value="For Team Meeting"/>	

### Event Description for Website Publication:

### Clinical Event Information:

Number of Incidents in Report:	Contact:	Alternative Person Title:	Alternative Person First Name:
<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Alternative Person Surname:	Alternative Person Phone:	Alternative Person Fax:	Alternative Person Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

### Recorded Problems Observed

#### Recorded Problems Observed:

### Clinical Signs, Symptoms and Conditions

#### Recorded Clinical Signs, Symptoms and Conditions:

### Health Impact

#### Recorded Health Impacts:

### Patient Information

Sex:	Weight:	Age:
<input type="text" value="S22"/>	<input type="text"/>	<input type="text" value="S22"/>
Patient Focused Corrective Action Taken:	Patient History:	
<input type="text"/>	<input type="text"/>	
Patient Outcome/Consequences:	Additional Event Description:	
<input type="text" value="Patient passed away."/>	<input type="text"/>	
Describe any test (Lab, xray, etc.):	Injured - Extent of Injury:	Other medical devices currently using/implanted:
<input type="text"/>	<input type="text"/>	<input type="text"/>
		Medical Problem Device Used For:
		<input type="text"/>

Add tional Patients Added:

0

Submitting Reporter Section

Search Reporter By Surname:

§22

Reporter Title:

§22

Position:

§22

Address 1:

2 Alma Road, Macquarie Park, NSW, 2113

Country:

Australia

Mobile:

Reporter #:

First Name:

§22

Address 2:

Postcode:

2113

Email:

§22 @medtronic.com

Surname:

§22

Company/Inst tution:

Medtron c Australasia Pty Ltd.

Town/Suburb:

Macquarie Park

Phone:

§22

Last External Submission By:

114395\_837 - 07/07/2022 14:50

Preferred Contact Method:

State:

NSW

Fax:

In tial Reporter Section

As Above?:

No

Search Reporter By Surname:

Title:

Position:

Address 1:

Postcode:

Mobile:

If No, fill out the following:

In tial Reporter #:

First Name:

Address 2:

Country:

Email:

Surname:

Company/Inst tution:

ST. ANDREWS

Town/Suburb:

Adelaide

Phone:

Allow the dev ce company to contact you about the incident:

Initial Reporter Confidential:

Yes

Preferred Contact Method:

State:

SA

Fax:

Dev ce Information Sect on

Product Exempt (Note: If not exempt, enter ARTG No):

No

Product Licence Category:

Included

Brand Name:

EEA - Intraluminal circular stapler, single-use

Model #:

TRIEEA28MT

Purchase Date:

Date of Inital Procedure:

Access Contact First Name:

Search Device ARTG:

293808

Device Class:

Class IIb

In tial Device Description:

EEA - Intraluminal circular stapler, single-use

Serial #:

Expiry Date:

Place of Implantation:

Access Contact Surname:

Device ARTG #:

293808

GMDN / UMDN Code:

59875

Usage of Device:

Batch #:

Date of Implant:

Reported Dev ce Locat on:

Access Contact Phone:

Therapeut c Licence Type:

Med cal Device

GMDN / UMDN Text:

Intraluminal circular stapler, single-use

Software Version:

Lot #:

P1L0118

Date of Explant:

Access Contact T tle:

Access Contact Fax:

<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Access Contact Email:	Licence Status:	Status Effective Date:	Additional Devices Added:
<input type="text"/>	A	13/09/2017	0

Manufacturer Information Section

Manufacturer Name:	Manufacturer Client Id:	Address 1:
Covidien Ilc	54968	<input type="text"/>
Address 2:	State/Province:	Country:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Postcode:	Fax:	Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Manufacturer Informed:	Contact Title:	Contact First Name:
Yes	<input type="text"/>	<input type="text"/>
Contact Surname:		
<input type="text"/>		

Supplier Information Section

Supplier Name:	Address 1:	Address 2:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Town/Suburb:	Country:	Postcode:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Phone:	Email:	Website:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Supplier Informed:	Contact Title:	Contact First Name:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Contact Surname:	Contact Fax:	Contact Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>

Report Status

For website publication:	Ready for Publication:	Investigated:	Investigation Reason:	Team Assignment:	Team Priority:
Yes	Yes	No	Rate considered low at this stage	Team B (Iib, Is & Im)	Routine

Team Review

Reviewed by Team:	Reason Sent To Meeting:	Outcome from team meeting:
<input type="text"/>	Death	Post-market Review is recommended
Notes for Team meeting:		
0 similar events, 0 PMR's, 0 recalls, sent to team meeting because of death. Should we ask whether or not the device was returned and why they haven't conducted any analysis if the device was returned?		
Outcomes from Team Meeting:		
Closed as this is a known procedural complication.		

In trial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	Assessor:
11/07/2022	9 - Death	<input type="text"/>	Yes	\$22
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:	Injured Party:
<input type="text"/>	No	Level 2 Investigation (for a single DIR)	Level 2 Investigation (for a single DIR)	Patient
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:	Found Prior To Use:	Manufacturer documentation:
<input type="text"/>				Unknown - updated information from the manufacturer is required
				Device Recalls:
				0. No recalls for similar incidents in Australia
				Is AE covered by current recall:



Investigat on Information Sect on - Submitted by Sponsor/Manufacturer

**Device Analysis Results:**

Medtronic conducted an investigat on based upon all received informat on.

The following complaint(s) were investigated:  
 It was reported that:  
 There was difficulty removing the dev ce from the patient.  
 - There was not enough information to make any determinat on based on insufficient evidence. The most likely root cause was not established.  
 The anvil disengaged either fully or partially from the device.  
 - There was not enough information to make any determinat on based on insufficient evidence. The most likely root cause was not established.

A medical assessment was performed. The medical assessment results were inconclusive.  
 A device history review was performed. The device history review determined there was no association to the complaint.

**Corrective/Preventative Actions:**

Further action was not required because there was insuff cient informat on.


**Details of Similar Events:**

**Australia:**  
 2 similar events from 3993 un ts (0.05009%) sold in CY2019;  
 2 similar events from 2977 un ts (0.06718%) sold in CY2020;  
 1 similar events from 3517 un ts (0.02843%) sold in CY2021;  
 1 similar events from 1580 un ts (0.06329%) sold in CY2022 (Includes this event).

**Worldw de:**  
 17 similar events from 46545 un ts (0.03652%) sold in CY2019;  
 9 similar events from 45149 units (0.01993%) sold in CY2020;  
 9 similar events from 71491 units (0.01259%) sold in CY2021;  
 9 similar events from 37992 units (0.02369%) sold in CY2022.

Similar inc dent rates were determined using FDD annex A code (DIFFICULT TO REMOVE - A150207) for the platform family type of product reported per subject event.

**Add tional Details (use for tables):**



**CAPA# Reference:**

**Risk Assessment**

<b>Frequency:</b> <input type="text"/>	<b>Sever ty:</b> <input type="text"/>
<b>Rating:</b> <input type="text"/>	<b>Type Cause and Outcome:</b> <input type="text"/>
<b>Expected Rate:</b> <input type="text"/>	<b>Number of Similar Events:</b> <input type="text"/>
<b>Actual Rate:</b> <input type="text"/>	

**Countries Similar Events Also Occurred:**

Austria, Finland, France, Germany, Ireland, Italy, Japan, Macao, Netherlands, New Zealand, Spain, Sweden, Switzerland, United Kingdom & United States

**Completed Actions:**

**Planned Act ons and Proposed Timelines:**

**Addit onal Comments:**

DIR Closed on 19 July 2022.

**Reason for Level 1 Investigation**

<b>Details of Reasons</b>	
<b>Reason for Level 1 Investigat on</b>	

**Focus of Level 2 Investigation**

<b>Details of Focus</b>	
<b>Essential Principles</b>	<b>If 'Other' Selected</b>

Sources of Evidence for Level 2

Details of Source			
Sources of Evidence	If 'Others' please specify here	Expected Sourcing Date	Date of Evidence Received

Evidence

Investigation Questions (Level 1 and Level 2):

Potential Risks

Delays in response by product manufacturers:

Delays in response by incident reporters:

Delays in analysis within the TGA:

Delays in reporting by other sources (e.g. clinical registries):

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:

Next Steps for Level 2 Investigation:

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Mechanical Problem	Retract on Problem		
Mechanical Problem	Detachment of Device or device Component		

Clinical signs symptoms and conditions

Details

Level 1	Level 2	Level 3
Infections	Sepsis	
General Disorders	Multiple Organ Dysfunction Syndrome	

Health Impact

Details

Level 1	Level 2	Level 3
Surgical Intervention	Additional Surgery	
Death		

Investigation Findings

Finding Details

Investigation Findings (Level 1)	Investigat on Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
No Findings Available			

Investigation Conclusion

Conclus on Details

Investigation Conclus on (L1)	Investigation Conclusion (L2)	If Add tional Conclus on Detail Requested
Cause Not Established		

Investigation Outcomes

Outcome Details

Outcome of Investigation (L1)	Outcome of Investigat on (L2)	If Add tional Conclus on Detail Requested
Reviewed, for Trending Purposes Only		

Investigation Summary

Latest Investigat on (DII) where this DIR is the Primary DIR: <input type="text"/>	Latest Investigat on (DII) where this DIR is a Related DIR: <input type="text"/>	Investigator: <input type="text"/>	Peer Review: <input type="text" value="No"/>	
Investigator's Notes: <input type="text"/>		Summary Findings: <input type="text" value="No further investigation will occur at this time; however the TGA will continue to mon tor the rate and pattern of occurrence and may re-open the file as appropriate."/>		Recall Number: <input type="text"/>

**Note:** Letter generat on buttons disabled if report not ready for webs te publ cation or risk analysis not completed.

Device Lookup

This section is used to match informat on prov ded via UDIR forms to ARTG information. You can select a Brand/Name from informat on prov ded in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered): <input type="text"/>	Brand Name: <input type="text"/>	Manufacturer Name: <input type="text"/>	Device ARTG #: <input type="text"/>
---	-------------------------------------	--	--

Other Devices

Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Cl ck **New** to begin entering informat on.

Rec No
1

Samples Record - Cl ck **[N]** to begin entering information. **Note:** Sample # Generated on Save.



Rec No	Details	Sample Details			Additional Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Device Description:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Version Number:	
					Who sent the device to the TGA?:			Why does the TGA have the sample?:	

Additional Patients

Click [N] to begin entering information.

Patient Details			
Sex:	Weight:	Age:	
Patient Focused Corrective Act on Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

Additional Device Information

Where did you get this device from?:  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

Similar Events

Similar events - how many times?:  Date of Recent Report:  Event Reported To:  Reporter Reference Number:

Device Access - Alternate Device Contact Information Provided

Title:  First Name:  Last Name:  Phone:

Fax:  Email:

Incident Location Details

Occurred in Australia:  Organisation:  Address Line 1:  Address Line 2:

Town/Suburb:  State:  Postcode:



## Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
335428	DIR-REQ		Closed	<b>S22</b>	OPR Administration User	19/07/2022	Normal	0

## Signature Details

Role	IRIS Investigator	
User	<b>S22</b>	
Signed At	19/07/2022 14:12:18	
Comment		



## Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 46 - ID : 547498

Released by **S22** on 28/09/2022 14:44:45

Report #:  Records Management #:  Reporter's Reference #:  Report Type:

ARTG: 293808

[Document Container URL](#)

### Report Information Section

Report Status:  Sponsor's Reported Category:  Date of Adverse Event:  Date of Initial Report:

Date of Final Report:  Date of Initial TGA Action:  Reviewed by Team:  Date Response Received:

Date Completed:  Operator at Time of Event:  If 'Other' Operator Selected:  Reporter consents to contact by sponsor:

Source of Report:  If 'Other' Source Selected:  Type of Initial Act on:

### Event Description for Website Publication:

### Clinical Event Information:

According to the reporter, two days post-operatively of laparoscopic low anterior bowel resection, the staple line leaked and had to have a second procedure. The patient is still in hospital and there is a high potential that the bowel cannot be re-anastomosed resulting in permanent colostomy bag for the patient. Further bowel needed to be resected due to the circular staple line dehiscence and it is unknown there will be enough bowel to rejoin. But during the first surgery, the device that was used was functioned as expected, no device misfired, air leak was performed successfully, and distal and proximal donuts were both intact. A second procedure was performed which is Hartmann's procedure and colostomy bag. These were to close rectal stump and create colostomy. Another operation is required but patient will undergo chemotherapy and have a brief healing period before further surgeries required.

Number of Incidents in Report:  Contact:  Alternative Person Title:  Alternative Person First Name:

Alternative Person Surname:  Alternative Person Phone:  Alternative Person Fax:  Alternative Person Email:

### Recorded Problems Observed

#### Recorded Problems Observed:

### Clinical Signs, Symptoms and Conditions

#### Recorded Clinical Signs, Symptoms and Conditions:

### Health Impact

#### Recorded Health Impacts:

### Patient Information

Sex:  Weight:  Age:

Patient Focused Corrective Action Taken:  Patient History:

Patient Outcome/Consequences:  Additional Event Description:

Describe any test (Lab, xray, etc.):

Injured - Extent of Injury:

Other medical devices currently using/implanted:

Medical Problem Device Used For:

Add tional Patients Added:

0

Submitting Reporter Section

Search Reporter By Surname:

s22

Reporter Title:

s22

Position:

s22

Address 1:

2 Alma Road

Country:

Australia

Mobile:

Reporter #:

First Name:

s22

Address 2:

Postcode:

2113

Email:

s22@medtron c.com

Surname:

s22

Company/Institut on:

Medtronic Australasia Pty Ltd

Town/Suburb:

Macquarie Park

Phone:

s22

Last External Submiss on By:

113291\_837 - 02/02/2023 12:45

Preferred Contact Method:

In tial Reporter Section

As Above?:

No

Search Reporter By Surname:

s47F

Title:

s47F

Position:

s47F

Address 1:

s47F

Postcode:

s47

Mobile:

If No, fill out the following:

In tial Reporter #:

First Name:

Address 2:

Country:

Email:

s47F

Surname:

Company/Institut on:

s47F

Town/Suburb:

s47F

Phone:

s47F

Allow the device company to contact you about the inc dent:

Initial Reporter Conf dential:

Yes

Preferred Contact Method:

Dev ce Information Sect on

Product Exempt (Note: If not exempt, enter ARTG No):

No

Product Licence Category:

Included

Brand Name:

EEA - Intraluminal circular stapler, single-use

Model #:

TRIEEA28XT

Purchase Date:

s47F

Date of Inital Procedure:

s47F

Access Contact First Name:

Search Device ARTG #:

293808

Dev ce Class:

Class IIb

In tial Device Descript on:

EEA - Intraluminal circular stapler, single-use

Serial #:

s47F

Expiry Date:

s47F

Place of Implantation:

s47F

Access Contact Surname:

Device ARTG #:

293808

GMDN / UMDN Code:

59875

Usage of Dev ce:

Single Use

Batch #:

s47F

Date of Implant:

s47F

Reported Device Location:

s47F

Access Contact Phone:

Therapeutic L cence Type:

Medical Device

GMDN / UMDN Text:

Intraluminal circular stapler, single-use

Software Vers on:

s47F

Lot #:

POC1305Y

Date of Explant:

s47F

Access Contact Title:

s47F

Access Contact Fax:

<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Access Contact Email:	Licence Status:	Status Effective Date:	Addit onal Devices Added:
<input type="text"/>	A	13/09/2017	0

Manufacturer Informat on Section

Manufacturer Name:	Manufacturer Client Id:	Address 1:
Covidien Ilc	54968	<input type="text"/>
Address 2:	State/Province:	Country:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Postcode:	Phone:	Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Manufacturer Informed:	Date Aware of Adverse Event:	Contact First Name:
Yes	12/10/2022	<input type="text"/>
Contact Surname:		
<input type="text"/>		

Supplier Informat on Section

Supplier Name:	Address 1:	Address 2:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Town/Suburb:	Country:	Postcode:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Phone:	Email:	Website:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Supplier Informed:	Contact Title:	Contact First Name:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Contact Surname:	Contact Fax:	Contact Email:
<input type="text"/>	<input type="text"/>	<input type="text"/>

Report Status

For website publ cat on:	Ready for Publication:	Investigated:	Investigation Reason:	Team Assignment:	Team Pr ority:
Yes	Yes	No	Device not returned	Team B (Iib, Is & Im)	Not Investigated

Team Review

Reviewed by Team:	Reason Sent To Meeting:	Outcome from team meeting:
<input type="text"/>	<input type="text"/>	<input type="text"/>
Notes for Team meeting:		
<input type="text"/>		
Outcomes from Team Meeting:		
<input type="text"/>		

In tial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
Date:	Sever ty:	Incidents in the last 12 months:	Manufacturer analysis:	Assessor:
09/02/2023	6 - Ongoing minor impairment	<input type="text"/>	Yes	\$22
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:	Injured Party:
<input type="text"/>	Unable to take action (device was not returned for analysis)	Level 1 Investigat on (to complete screening)	Screening only	Patient
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:	Found Prior To Use:	Manufacturer documentat on:
				Unknown - updated informat on from the manufacturer is required
				Device Recalls:
				0. No recalls for similar inc dents in Australia
				Is AE covered by current recall:

<input type="text"/>	<input type="text"/>	The dev ce was not returned - Low rate		No	<b>Document 6</b>
Incidents Worldw de: <input type="text"/>	Number of potential contributing factors: No			Reusable: No	Similar events (past 6 months): 0 incidents
Products supplied the last 12 months: <input type="text"/>	Specific factors identified: <input type="text"/>	ESTIMATED LEVEL OF PRIORITY: Routine	FINAL LEVEL OF PRIORITY: Routine		
Products supplied last 24 months: <input type="text"/>	Number of potential sens tivities: No	EXCEPTION TO PRIORITY LEVEL: <input type="text"/>		3 or more events - batch/model: <input type="text"/>	
Products supplied last 36 months: <input type="text"/>	Specific sens tivities identified: <input type="text"/>			3 or more events - health district: <input type="text"/>	
Products supplied Worldw de: <input type="text"/>	Consultations during risk assessment: I did none of the above incidents	Final Risk Assessment: Yes	3 or more events - organisation: <input type="text"/>		

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section				
Update Device Details?:	<b>Background Information</b>	<b>Risk Assessment - Section A</b>	<b>Risk Assessment - Section B</b>	<b>Risk Assessment - Section C</b>	<b>Risk Assessment - Section D</b>
	Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	
Copy Data From:					
	Assessor:	Manufacturer documentation:	Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION: FINAL LEVEL OF INVESTIGATION:
	Injured Party:	Device Recalls:	Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:
	Found Prior To Use:	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors:	
	Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY: FINAL LEVEL OF PRIORITY:
		3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:
		3 or more events - health district:	Products supplied last 36 months:	Specific sensitivities identified:	
		3 or more events - organisation:	Products supplied Worldw de:	Consultations during risk assessment:	Final Risk Assessment:

Sponsor/Manufacturer Information Section

Search Sponsors: 837	Name: Medtronic Australasia Pty Ltd	Client #: 837
Attention To: <input type="text"/>	Address 1: PO Box 945	Address 2: <input type="text"/>
State: NSW	Postcode: 1670	Phone: <input type="text"/>
Email: <input type="text"/>		Town/Suburb: NORTH RYDE BC
		Fax: <input type="text"/>

Investigat on Information Sect on - Submitted by Sponsor/Manufacturer

Device Analysis Results:

Medtronic conducted an investigat on based upon all received informat on.

The following complaint(s) were investigated:

It was reported that:

An unexpected content leak occurred along the staple line.

- There was not enough information to make any determinat on based on insufficient evidence. The most likely root cause was not established.

An incision of >1 inch was needed resulting from product failure.

- There was not enough information to make any determinat on based on insufficient evidence. The most likely root cause was not established.

A device history review was performed. The device history review determined there was no association to the complaint.

Further action was not required because there was insufficient information.

Should new informat on become available the file will be re-opened and the investigation summary will be amended as appropriate.

Addit onal Details (use for tables):



Details of Similar Events:

Australia:

0 similar events from 379 un ts (0%) sold in CY2020;

1 similar events from 567 un ts (0.17637%) sold in CY2021;

3 similar events from 294 un ts (1.02041%) sold in CY2022;

0 similar events from 81 units (0%) sold in CY2023 (Includes this event).

Worldw de:

14 similar events from 6400 un ts (0.21875%) sold in CY2020;

32 similar events from 14244 units (0.22466%) sold in CY2021;

64 similar events from 17296 units (0.37003%) sold in CY2022;

1 similar events from 2535 units (0.03945%) sold in CY2023.

Similar inc dent rates were determined using FDD annex A code (ADVERSE EVENT WITHOUT IDENTIFIED DEVICE OR USE PROBLEM - A24) for the platform family type of product reported per subject event.

CAPA# Reference:

Risk Assessment

Frequency:

Severity:

Rating:

Type Cause and Outcome:

Number of Similar Events:

Expected Rate:

Actual Rate:

Countries Similar Events Also Occurred:

Australia, Belgium, Brazil, Canada, Chile, Finland, France, Germany, Greece, Hong Kong, India, Ireland, Israel, Italy, Japan, Korea Republic Of, Malaysia, Mexico, Netherlands, Norway, Poland, South Afr ca, Spain, Sw tzerland, Thailand, Un ted Kingdom, Un ted States

Completed Actions:

Planned Actions and Proposed Timelines:

Addit onal Comments:

Similar inc dent rates were determined using FDD annex A code (DEVICE ALARM SYSTEM - A1601) for the platform family type of product reported per subject event.

DIR closed on 9/2/2023

Reason for Level 1 Investigation

Details of Reasons

Reason for Level 1 Investigat on

Focus of Level 2 Investigation

Details of Focus

Essential Principles

If 'Other' Selected

Sources of Evidence for Level 2

Details of Source

Sources of Ev dence

If 'Others' please specify here

Expected Sourcing Date

Date of Ev dence Received



Evidence

Investigation Questions (Level 1 and Level 2):

Potential Risks

Delays in response by product manufacturers:

Delays in response by incident reporters:

Delays in analysis within the TGA:

Delays in reporting by other sources (e.g. clinical registries):

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:

Next Steps for Level 2 Investigation:

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Mechanical Problem	Leak / Splash	Fluid Leak	

Clinical signs symptoms and conditions

Details

Level 1	Level 2	Level 3
Gastrointestinal System	Rectal Anastomotic Leakage	

Health Impact

Details

Level 1	Level 2	Level 3
Surgical Intervention	Additional Surgery	

Investigation Findings

Finding Details

Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
No Findings Available			

Investigation Conclusion

Conclus on Details

Investigation Conclus on (L1)	Investigation Conclusion (L2)	If Add tional Conclus on Detail Requested
Cause Not Established		

Investigation Outcomes

Outcome Details		
Outcome of Investigation (L1)	Outcome of Investigat on (L2)	If Add tional Conclus on Detail Requested
Reviewed, for Trending Purposes Only		

Investigation Summary

Latest Investigat on (DII) where this DIR is the Primary DIR: <input type="text"/>	Latest Investigat on (DII) where this DIR is a Related DIR: <input type="text"/>	Investigator: <input type="text"/>	Peer Review: <input type="text" value="No"/>	
Investigator's Notes: <input type="text"/>		Summary Findings: <input type="text" value="No further investigation will occur at this time; however, the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate."/>		Recall Number: <input type="text"/>

**Note:** Letter generat on buttons disabled if report not ready for webs te publ cation or risk analysis not completed.

Device Lookup

This section is used to match informat on prov ded via UDIR forms to ARTG information. You can select a Brand/Name from informat on prov ded in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered): <input type="text"/>	Brand Name: <input type="text"/>	Manufacturer Name: <input type="text"/>	Device ARTG #: <input type="text"/>
---	-------------------------------------	--	--

Other Devices

Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Cl ck **New** to begin entering informat on.

Rec No
1

Samples Record - Cl ck **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details			Add t onal Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Dev ce Descript on:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Vers on Number:	
						Who sent the dev ce to the TGA?:		Why does the TGA have the sample?:	

Additional Patients

Click [N] to begin entering information.

Patient Details			
Sex:	Weight:	Age:	
Patient Focused Corrective Act on Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

Additional Device Information

Where did you get this device from?:  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

Similar Events

Similar events - how many times?:  Date of Recent Report:  Event Reported To:  Reporter Reference Number:

Device Access - Alternate Device Contact Information Provided

Title:  First Name:  Last Name:  Phone:


Fax:  Email:

Incident Location Details

Occurred in Australia:  Organisation:  Address Line 1:  Address Line 2:

Town/Suburb:  State:  Postcode:

Flow Details DIR-REQ - Device Incident Request 352894

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
352894	DIR-REQ		Closed		OPR Administration User	09/02/2023	Normal	0

Signature Details

Role	IRIS Investigator	
User	522	
Signed At	09/02/2023 15:47:26	
Comment		



# Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 46 - ID : 564527

Released by **S22** on 28/09/2022 14:44:45

Report #: <input type="text" value="87436"/>	Records Management #: <input type="text"/>	Reporter's Reference #: <input type="text" value="705600221"/>	Report Type: <input type="text" value="Final"/>
---	---	---	--

ARTG: 293808	<a href="#">Document Container URL</a>		
--------------	--	--	--

Report Information Section

Report Status: <input type="text" value="Closed"/>	Sponsor's Reported Category: <input type="text"/>	Date of Adverse Event: <input type="text" value="S22"/>	Date of Initial Report: <input type="text" value="10/05/2023"/>
Date of Final Report: <input type="text" value="04/07/2023"/>	Date of Initial TGA Action: <input type="text" value="10/05/2023"/>	Reviewed by Team: <input type="text"/>	Date Response Received: <input type="text"/>
Date Completed: <input type="text" value="05/07/2023"/>	Operator at Time of Event: <input type="text"/>	If 'Other' Operator Selected: <input type="text"/>	Reporter consents to contact by sponsor: <input type="text" value="N/A"/>
Source of Report: <input type="text" value="Sponsor"/>	If 'Other' Source Selected: <input type="text"/>	Type of Initial Action: <input type="text" value="Trend data only"/>	

Event Description for Website Publication:

Clinical Event Information:

Number of Incidents in Report: <input type="text" value="1"/>	Contact: <input type="text"/>	Alternative Person Title: <input type="text"/>	Alternative Person First Name: <input type="text"/>
Alternative Person Surname: <input type="text"/>	Alternative Person Phone: <input type="text"/>	Alternative Person Fax: <input type="text"/>	Alternative Person Email: <input type="text"/>

Recorded Problems Observed

Recorded Problems Observed:

Clinical Signs, Symptoms and Conditions

Recorded Clinical Signs, Symptoms and Conditions:

Health Impact

Recorded Health Impacts:

Patient Information

Sex: <input type="text"/>	Weight: <input type="text"/>	Age: <input type="text"/>
Patient Focused Corrective Action Taken: <input type="text"/>	Patient History: <input type="text"/>	
Patient Outcome/Consequences: <input type="text"/>	Additional Event Description: <input type="text"/>	
Describe any test (Lab, xray, etc.): <input type="text"/>	Injured - Extent of Injury: <input type="text"/>	Other medical devices currently using/implanted: <input type="text"/>
		Medical Problem Device Used For: <input type="text"/>

Additional Patients Added:

ADD OTHER PATIENTS ABOVE.

0

Submitting Reporter Section

Search Reporter By Surname:

Reporter Title:

Position:

Address 1:

Country:

Mobile:

Reporter #:

First Name:

Address 2:

Postcode:

Email:

Surname:

Company/Inst tution:

Town/Suburb:

Phone:

Last External Submission By:

Preferred Contact Method:

State:

Fax:

In tial Reporter Section

As Above?:

If No, fill out the following:

Initial Reporter Confidential:

Search Reporter By Surname:

In tial Reporter #:

Title:

First Name:

Surname:

Position:

Company/Inst tution:

Address 1:

Address 2:

Town/Suburb:

State:

Postcode:

Country:

Phone:

Fax:

Mobile:

Email:

Allow the dev ce company to contact you about the inc dent:

Dev ce Information Sect on

Product Exempt (Note: If not exempt, enter ARTG No):

Search Device ARTG:

Device ARTG #:

Therapeut c Licence Type:

Product Licence Category:

Device Class:

GMDN / UMDN Code:

GMDN / UMDN Text:

Brand Name:

In tial Device Description:

Usage of Device:

Software Version:

Model #:

Serial #:

Batch #:

Lot #:

Purchase Date:

Expiry Date:

Date of Implant:

Date of Explant:

Date of Inital Procedure:

Place of Implantat on:

Reported Dev ce Locat on:

Access Contact T tle:

Access Contact First Name:

Access Contact Surname:

Access Contact Phone:

Access Contact Fax:

Access Contact Email: <input type="text"/>	Licence Status: <input type="text" value="A"/>	Status Effective Date: <input type="text" value="13/09/2017"/>	Add t onal Devices Added: <input type="text" value="0"/>
---	---	---	---

Manufacturer Informat on Section

Manufacturer Name: <input type="text" value="Covidien llc"/>	Manufacturer Client Id: <input type="text" value="54968"/>	Address 1: <input type="text"/>
Address 2: <input type="text"/>	State/Province: <input type="text"/>	Country: <input type="text"/>
Postcode: <input type="text"/>	Phone: <input type="text"/>	Fax: <input type="text"/>
Manufacturer Informed: <input type="text" value="Yes"/>	Date Aware of Adverse Event: <input type="text" value="05/05/2023"/>	Contact Title: <input type="text"/>
Contact Surname: <input type="text"/>		Contact First Name: <input type="text"/>

Supplier Informat on Section

Supplier Name: <input type="text"/>	Address 1: <input type="text"/>	Address 2: <input type="text"/>
Town/Suburb: <input type="text"/>	Country: <input type="text"/>	Postcode: <input type="text"/>
Phone: <input type="text"/>	Email: <input type="text"/>	Website: <input type="text"/>
Supplier Informed: <input type="text"/>	Contact Title: <input type="text"/>	Contact First Name: <input type="text"/>
Contact Surname: <input type="text"/>	Contact Fax: <input type="text"/>	Contact Email: <input type="text"/>

Report Status					
For website publ cat on: <input type="text" value="Yes"/>	Ready for Publication: <input type="text" value="Yes"/>	Investigated: <input type="text" value="No"/>	Investigation Reason: <input type="text" value="Event determined to be most likely due to user error"/>	Team Assignment: <input type="text" value="Team B (Iib, Is &amp; Im)"/>	Team Pr ority: <input type="text" value="Not Investigated"/>

Team Review			
Reviewed by Team: <input type="text"/>	Reason Sent To Meeting: <input type="text"/>	Outcome from team meeting: <input type="text"/>	
Notes for Team meeting: <input type="text"/>			
Outcomes from Team Meeting: <input type="text"/>			

In tial Risk Analysis					
Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
Date: <input type="text" value="10/05/2023"/>	Sever ty: <input type="text" value="5 - An illness/injury was resolved or prevented with treatment by a health profess onal"/>	Incidents in the last 12 months: <input type="text"/>	Manufacturer analysis: <input type="text" value="Yes"/>	Assessor: <input type="text" value="S22"/>	Manufacturer documentat on: <input type="text" value="Unknown - updated informat on from the manufacturer is required"/>
Incidents in last 24 months: <input type="text"/>	Manufacturer action: <input type="text" value="No"/>	ESTIMATED LEVEL OF INVESTIGATION: <input type="text" value="Level 1 Investigat on (to complete screening)"/>	FINAL LEVEL OF INVESTIGATION: <input type="text" value="Screening only"/>	Injured Party: <input type="text" value="Patient"/>	Device Recalls: <input type="text" value="0. No recalls for similar inc dents in Australia"/>
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:	Found Prior To Use:	Is AE covered by current recall:	

<input type="text"/>	<input type="text"/>	<input type="text" value="The most likely root cause was traced to the user."/>		<input type="text" value="No"/>	<input type="text" value="Document 7"/>
Incidents Worldw de: <input type="text"/>	Number of potential contributing factors: <input type="text" value="No"/>			Reusable: <input type="text" value="No"/>	Similar events (past 6 months): <input type="text" value="0 incidents"/>
Products supplied the last 12 months: <input type="text"/>	Specific factors identified: <input type="text"/>	ESTIMATED LEVEL OF PRIORITY: <input type="text" value="Routine"/>	FINAL LEVEL OF PRIORITY: <input type="text" value="Routine"/>	3 or more events - batch/model: <input type="text" value="No"/>	
Products supplied last 24 months: <input type="text"/>	Number of potential sens tivities: <input type="text" value="No"/>	EXCEPTION TO PRIORITY LEVEL: <input type="text" value="The most likely root cause was traced to the user."/>		3 or more events - health district: <input type="text" value="No"/>	
Products supplied last 36 months: <input type="text"/>	Specific sens tivities identified: <input type="text"/>			3 or more events - organisation: <input type="text" value="No"/>	
Products supplied Worldw de: <input type="text"/>	Consultations during risk assessment: <input type="text" value="I did none of the above incidents"/>	Final Risk Assessment: <input type="text" value="Yes"/>			

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section				
Update Device Details?:	<b>Background Information</b>	<b>Risk Assessment - Section A</b>	<b>Risk Assessment - Section B</b>	<b>Risk Assessment - Section C</b>	<b>Risk Assessment - Section D</b>
	Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	
Copy Data From:					
	Assessor:	Manufacturer documentation:	Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION: FINAL LEVEL OF INVESTIGATION:
	Injured Party:	Device Recalls:	Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:
	Found Prior To Use:	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors:	
	Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY: FINAL LEVEL OF PRIORITY:
		3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:
		3 or more events - health district:	Products supplied last 36 months:	Specific sensitivities identified:	
		3 or more events - organisation:	Products supplied Worldw de:	Consultations during risk assessment:	Final Risk Assessment:

Sponsor/Manufacturer Information Section

Search Sponsors: <input type="text" value="837"/>	Name: <input type="text" value="Medtronic Australasia Pty Ltd"/>	Client #: <input type="text" value="837"/>
Attention To: <input type="text"/>	Address 1: <input type="text" value="PO Box 945"/>	Address 2: <input type="text"/>
State: <input type="text" value="NSW"/>	Postcode: <input type="text" value="1670"/>	Town/Suburb: <input type="text" value="NORTH RYDE BC"/>
Email: <input type="text"/>	Phone: <input type="text"/>	Fax: <input type="text"/>



Investigat on Information Sect on - Submitted by Sponsor/Manufacturer

Device Analysis Results:

Medtronic conducted an investigat on based upon all received informat on.

The following complaint(s) were investigated:  
 It was reported that:  
 - The staples d d not deploy from the stapler after the handle was squeezed.  
 - This complaint was confirmed based on the actual device. The most likely root cause was traced to the user.

A medical assessment was performed. The medical assessment determined that there was a relat on to the complaint.

Product analysis occurred.  
 The primary analysis finding was:  
 - The anvil was tilted. Microscopic inspect on of the anvil mylar ring noted that t was crushed.

Further action was not required because the event had foreseen risk and is included in a data monitoring plan.

Should new informat on become available the file will be re-opened and the investigation summary will be amended as appropriate.

Corrective/Preventative Actions:

Further action was not required because the event had foreseen risk and is included in a data mon toring plan.  
 Should new informat on become available the file will be re-opened and the investigat on summary will be amended as appropriate.

Details of Similar Events:

Australia:  
 #DIV/0!  
 0 similar events from 124 units (0%) sold in CY2021;  
 0 similar events from 257 units (0%) sold in CY2022;  
 1 similar events from 109 units (0.91743%) sold in CY2023 (Includes this event).

Worldw de:  
 1 similar events from 81 un ts (1.23457%) sold in CY2020;  
 2 similar events from 855 units (0.23392%) sold in CY2021;  
 7 similar events from 1622 un ts (0.43157%) sold in CY2022;  
 6 similar events from 1271 un ts (0.47207%) sold in CY2023.

Similar inc dent rates were determined using FDD annex A code (MECHANICS ALTERED - A0507) for the platform family type of product reported per subject event.

Add tional Details (use for tables):



CAPA# Reference:

Risk Assessment

Frequency:	Sever ty:
<input type="text"/>	<input type="text"/>
Rating:	
<input type="text"/>	
Expected Rate:	Actual Rate:
<input type="text"/>	<input type="text"/>

Type Cause and Outcome:	Number of Similar Events:
<input type="text"/>	<input type="text"/>

Countries Similar Events Also Occurred:

Australia, Belgium, Brazil, France, Germany, Italy, Japan, Saudi Arabia, Singapore, Switzerland, Thailand, United States

Completed Actions:

Planned Act ons and Proposed Timelines:

Addit onal Comments:

DIR closed on 05/07/2023.

Reason for Level 1 Investigation

Details of Reasons	
Reason for Level 1 Investigat on	

Focus of Level 2 Investigation

Details of Focus	
Essential Principles	If 'Other' Selected

Sources of Evidence for Level 2

Details of Source			
Sources of Evidence	If 'Others' please specify here	Expected Sourcing Date	Date of Evidence Received

Evidence

Investigation Questions (Level 1 and Level 2):

Potential Risks

Delays in response by product manufacturers:      
 Delays in response by incident reporters:      
 Delays in analysis within the TGA:      
 Delays in reporting by other sources (e.g. clinical registries):

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:      
 Next Steps for Level 2 Investigation:

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details									
Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details			
Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Mechanical Problem	Firing Problem	Failure to Fire	

Clinical signs symptoms and conditions

Details		
Level 1	Level 2	Level 3
Procedural Complications	Appropriate Term / Code Not Available	

Health Impact

Details		
Level 1	Level 2	Level 3
Surgical Intervention	Prolonged surgery	

Investigation Findings

Finding Details			
Investigation Findings (Level 1)	Investigat on Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
Usage Problem Identified			Root cause was likely traced to the user

Investigation Conclusion

Conclus on Details		
Investigation Conclus on (L1)	Investigation Conclusion (L2)	If Add tional Conclus on Detail Requested
Cause Traced to User		

Investigation Outcomes

Outcome Details		
Outcome of Investigation (L1)	Outcome of Investigat on (L2)	If Add tional Conclus on Detail Requested
Reviewed, for Trending Purposes Only		

Investigation Summary

Latest Investigat on (DII) where this DIR is the Primary DIR: <input type="text"/>	Latest Investigat on (DII) where this DIR is a Related DIR: <input type="text"/>	Investigator: <input type="text"/>	Peer Review: <input type="text" value="No"/>	
Investigator's Notes: <input type="text"/>		Summary Findings: <input type="text" value="No further investigation will occur at this time; however, the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate."/>		Recall Number: <input type="text"/>

**Note:** Letter generat on buttons disabled if report not ready for webs te publ cation or risk analysis not completed.

Device Lookup

This section is used to match informat on prov ded via UDIR forms to ARTG information. You can select a Brand/Name from informat on prov ded in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered): <input type="text"/>	Brand Name: <input type="text"/>	Manufacturer Name: <input type="text"/>	Device ARTG #: <input type="text"/>
---	-------------------------------------	--	--

Other Devices					
Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Cl ck **New** to begin entering informat on.

Rec No
1

Samples Record - Click **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details	Add t onal Details

1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Device Description:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Version Number:	
	Who sent the device to the TGA?:							Why does the TGA have the sample?:	

Additional Patients

Click [N] to begin entering information.

Patient Details			
Sex:	Weight:	Age:	
Patient Focused Corrective Act on Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:
Other Consequence:	Describe any test (Lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For:

Additional Device Information

Where did you get this device from?  How reliant is the affected person on correct/safe operation of this device?:

Any other relevant information to aid assessing/investigating the incident?:

Similar Events

Similar events - how many times?: <input type="text"/>	Date of Recent Report: <input type="text"/>	Event Reported To: <input type="text"/>	Reporter Reference Number: <input type="text"/>
--	---	---	---

Device Access - Alternate Device Contact Information Provided

Title: <input type="text"/>	First Name: <input type="text"/>	Last Name: <input type="text"/>	Phone: <input type="text"/>
Fax: <input type="text"/>	Email: <input type="text"/>		

Incident Location Details

Occurred in Australia: <input type="text"/>	Organisation: <input type="text"/>	Address Line 1: <input type="text"/>	Address Line 2: <input type="text"/>
Town/Suburb: <input type="text"/>	State: <input type="text"/>	Postcode: <input type="text"/>	

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
383149	DIR-REQ		Closed	S22	OPR Administration User	05/07/2023	Normal	0

Signature Details

Role	IRIS Investigator	
User	S22	
Signed At	05/07/2023 16:06:44	
Comment		



# Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

06/05/2014

SIGNED

DIR : 20 - ID : 275697

Released by s22 on 25/06/2015 15:11:06

## Report Information Section

Report #: 34124	Records Management #: 	Reporter's Reference #: US201404-0480	Report Type: Final
Report Status: Closed	Sponsor's Reported Category: Other	Date of Adverse Event: s22	Date of Initial Report: 06/05/2014
Date of Final Report: 31/03/2015	Date of Initial TGA Action: 06/05/2014	Reviewed by DIRE: 	Date Response Received: 
Date Completed: 31/03/2015	Operator at Time of Event: Healthcare Professional	If 'Other' Operator Selected: 	Reporter Confidentiality: No
Source of Report: Sponsor	If 'Other' Source Selected: 	Type of Initial Action: Trend data only	

## Event Description for Website Publication:

On day 6 post-operatively, the patient developed peritonitis due to a suspected anastomotic leak after the first attempt at opening their bowels. The surgeon had to re-operate on the patient and noted that on one of the lateral sides of the anastomosis there was a leak (detected only by the air-leak test).

## Clinical Event Information:

The surgeon provided feedback regarding two (2) recent colectomies. One case was a laparoscopic anterior resection of rectum with reconstruction via end-to-end anastomosis using a Covidien size 28 stapler. The operation went well and the anastomosis leak tests using dye and air demonstrated no leak. The anastomotic donuts were noted to be complete. As the procedure was done entirely laparoscopically, no reinforcement of the staple line with sutures occurred. On day 6 post-operatively, the patient developed peritonitis due to a suspected anastomotic leak after the first attempt at opening her bowels. The surgeon had to re-operate on the patient and noted that on one of the lateral sides of the anastomosis there was a leak (detected only by the air-leak test). The surgeon noted that the bowel was well vascularized and there was no tissue tension. The patient required an ileostomy and prolonged rehabilitation. The surgeon has questioned whether the leak was device-related.

Number of Incidents in Report: 1	Contact: 	Alternative Person Title: 	Alternative Person First Name: 
Alternative Person Surname: 	Alternative Person Phone: 	Alternative Person Fax: 	

## Patient Information

Sex: 	Weight: 	Age: 	
----------	-------------	----------	--

## Patient Focused Corrective Action Taken:

--

Patient History:

Patient Outcome/Consequences:

Other Devices Involved:

## Submitting Reporter Section

Search Reporter By Surname:

Reporter #:

Reporter Title:

First Name:

Surname:

Position:

Company/Institution:

Address 1:

Address 2:

Town/Suburb:

State:

Country:

Postcode:

Phone:

Fax:

Mobile:

Email:

Last External Submission By:

## Initial Reporter Section

As Above?:

*If No, fill out the following:*

Initial Reporter Confidential:

Search Reporter By Surname:

Initial Reporter #:

Title:

First Name:

Surname:

Position:

Company/Institution:

Address 1:

Address 2:

Town/Suburb:

State:

Postcode:

Phone:

Fax:

Mobile:

Email:

## Device Information Section

Product Exempt:

No

*If No, fill out ARTG No:*

Search Device ARTG:

178517

Device ARTG #:

178517

Therapeutic Licence Type:

Medical Device

Product Licence Category:

Included

Device Class:

Class IIb

GMDN / UMDN Code:

45183

GMDN / UMDN Text:

Applier, surgical staple, cutting

Brand Name:

DST EEA 28MM Single-Use Stapler

Initial Device Description:

DST EEA 28MM Single-Use Stapler

Usage of Device:

Single Use

Software Version:

Model #:

EEA28

Serial #:

Batch #:

Lot #:

unknown

Purchase Date:

Expiry Date:

Date of Implant:

Date of Explant:

Reported Device Location:

Access Contact Title:

Access Contact First Name:

Access Contact Surname:

Access Contact Phone:

Access Contact Fax:

## Manufacturer Information Section

Manufacturer Name:

Covidien llc

Manufacturer Client Id:

54968

Address 1:

15 Hampshire Street

Address 2:

Town/Suburb:

Mansfield

State/Province:

MA

Country:

United States

Postcode:

02048

Phone:

§22

Fax:

Email:

australia-regulatory-affairs@covidien.com

Manufacturer Informed:

Yes

Date Aware of Adverse Event:

08/04/2014

Contact Title:

§22

Contact First Name:

§22

Contact Surname:

§22



Supplier Information Section			
Supplier Name:		Address 1:	Address 2:
<input type="text"/>		<input type="text"/>	<input type="text"/>
Town/Suburb:	State:	Postcode:	Phone:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Fax:	Email:		Supplier Informed:
<input type="text"/>	<input type="text"/>		<input type="text"/>
Date of Supplier Contact:	Contact Title:	Contact First Name:	Contact Surname:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Contact Phone:	Contact Fax:		
<input type="text"/>	<input type="text"/>		

Statistics Checklist Section				
Date:	Assessed By:	For website publication:	Ready for Publication:	Exclude report from DIRE:
<input type="text" value="31/03/2015"/>	<input type="text" value="s22"/>	<input type="text" value="Yes"/>	<input type="text" value="Yes"/>	<input checked="" type="checkbox"/>
Sample Received:	Sterile:	Reusable:	Single Use:	Potential Effect:
<input type="text" value="No"/>	<input type="text" value="Yes"/>	<input type="text" value="No"/>	<input type="text" value="Yes"/>	<input type="text" value="Serious Injury"/>
Actual Effect:	Injured Party:			Risk Frequency:
<input type="text" value="Serious Injury"/>	<input type="text" value="Patient"/>			<input type="text" value="Unlikely"/>
Risk Severity:	Risk Detectability:	Classification:	Investigated:	Date of DIRE Meeting:
<input type="text" value="Serious"/>	<input type="text" value="Unlikely"/>	<input type="text" value="Not Investigated"/>	<input type="text" value="No"/>	<input type="text"/>
DIRE Meeting Notes:				
<input type="text" value="This report has been reassessed due to age of the report, as the final report was not submitted by the sponsor. A search of the IRIS database has not shown an increasing trend in this type of event."/>				

Sponsor Information Section			
Search Sponsors:	Name:	Client #:	
<input type="text" value="Covidien"/>	<input type="text" value="Covidien Pty Ltd"/>	<input type="text" value="283"/>	
Attention To:	Address 1:	Address 2:	Town/Suburb:
<input type="text" value="s22"/>	<input type="text" value="Locked Bag 2020"/>	<input type="text"/>	<input type="text" value="LANE COVE DC"/>
State:	Postcode:	Phone:	Fax:
<input type="text" value="NSW"/>	<input type="text" value="2066"/>	<input type="text" value="s22"/>	<input type="text" value="s22"/>

Email:

## Investigation Information Section

Device Analysis Results:

Corrective/Preventative Actions:

Details of Similar Events:

Number of Similar Events:

Rate of Similar Events:

Countries Similar Events Also Occurred:

Additional Comments:

## Device Lookup

This section is used to match information provided via UDIR forms to ARTG information. You can select a Brand/Name from information provided in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered):

Brand Name:

Manufacturer Name:

Device ARTG #:

## Other Devices

Device ARTG No	Manufacturer Name	Sponsor/Supplier	Trade/Brand Name	Serial #	Model Number	GMDN / UMDN Text

Related DIR Information - Click **New** to begin entering information.

## Incident Details

DIR #	Brand Name	Reporter First Name	Reporter Surname	Company/Institution

Samples Record - Click **[N]** to begin entering information.

#### Sample Details

Sample #	Sample Requested	Sample Received	# Samples from Reporter	# Samples from Sponsor	Outcome of TGA's Testing

Correspondence Information Section. Note that the Correspondence Recipient will receive a notification if the Date Received is not filled in by the Date Expected.

#### Correspondence Details

Correspondence Type	Correspondence recipient (email)	Date Sent	Date Response Expected	Date Received	Sponsor's Response	Investigator's Notes

List of Problem Type Codes - Click **[N]** to begin entering information.

#### Type Details

Type of Problem (Level 1)	Type of Problem (Level 2)	If 'Other' Type Selected
Usability	Use of Device Issue	

#### Investigation Problem Causes

##### Cause Details

Cause of Problem (Level 1)	Cause of Problem (Level 2)	If 'Other' Cause Selected
Unable to confirm complaint	Device not returned	

#### Investigation Outcomes

##### Outcome Details

Outcome of Investigation	If Additional Outcome Detail Requested
Reviewed, for Trending Purposes Only	

Recall Number:

Investigation Summary:

No further investigation will occur at this time, however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.

**Flow Details : DIR-REQ - Device Incident Request : 46126****Request Details**

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
46126	DIR-REQ		Closed	s22	OPR Administration User	31/03/2015	Normal	0

**Signature Details**

Role	IRIS Investigator	
User	s22	
Signed At	31/03/2015 15:16:42	
Comment		