

| Equipment Type (NATA General Equipment Description | | Calibrated | | Verified | | Comment | LIMS Number/s |
|--|---|---------------|----------|---------------|-----------|--|---------------|
| | | Internal/Exte | Interval | Internal/Exte | Interval3 | | |
| HUMICAP humidity probe | HMP76 Humidity and temperature probe - Vaisala | External | 1 year | Internal | 6 months | This is a 2 in 1 T/RH sensor. Only the temperature is verified every 6 months. | 32223 |
| HUMICAP humidity probe | HMT331 Humidity and temperature transmitter - Vaisala | External | 1 year | Internal | 6 months | This is a 2 in 1 T/RH sensor. Only the temperature is verified every 6 months. | 32222 |
| Glass Microfibre Filter Media | Grade HE-1071 (a.k.a. green-line media) supplied by ATI | | | Internal | daily | used to verify particle size distribution and ATI 100Xs system performance. Reference charts | 10059 |
| Temperature | Fluke Platinum RTD | External | 1 year | Internal | 6 months | Used to 6 monthly check Testo H/T sensors | 32924, 32186 |
| Temperature | Fluke 1586a Temperature Logger | External | 2 years | Internal | 6 Months | Used to perform spatial on Binder enclosures. This used in conjunction with spatial probes 32855-32861 | 32848 |
| Temperature | PRT - probes | External | 2 years | Internal | 6 Months | Used to perform spatial on Binder enclosures | 32855-32861 |
| Set of weights | Masscal 2 kg - 5 g stainless steel | External | 3 years | | | Used to perform monthly and 6 monthly user checks of balances throughout labs | 5815 |