

9. Finished product expiry specification:

<u>Test</u>	<u>Specification</u>	<u>Reference</u>
Description	A clear, colourless or almost colourless solution, free from visible particulate.	Visual
Identification	[REDACTED]	[REDACTED]
Colour & Clarity of Solution	[REDACTED]	[REDACTED]
Assay	[REDACTED]	[REDACTED]
Related Substances	[REDACTED]	[REDACTED]
pH	[REDACTED]	[REDACTED]
Removable Volume	Not less than nominal fill	BP
Sterility	Complies with test	TGAL
Pyrogens	Complies with test	BP

10. Test methodology:

No special equipment required.

Chemical test provided, see Appendix 1.